

Internal Hexagon
Product Catalog
2025

Internal Hexagon
Product Catalog

NOT FOR USE IN THE USA AND CANADA

Products may not be available in some regulated markets. Please contact your local Dentsply Sirona sales office for current product assortment and availability.

TABLE OF CONTENTS



IMPLANT SYSTEMS	4 - 27.
SEVEN XD	6
M4	18



TOOLS AND KITS	28-63.
Step Pilot Drills	30
SEVEN Twist Drills	32
M4 Twist Drills	33
SEVEN and M4 Short Twist Drills	34
Marking Drills and Countersinks	36
Trephine Burrs	38
Insertion Tools	39
Surgical Tools	40
Prosthetic Instruments	42
SEVEN XD Placement Set	45
SEVEN Surgical Kit	47
SEVEN Drill Replacement Kit	51
M4 Surgical Kit	52
Bone Profiler Kit	55
SEVEN Drill Stoppers Kit, Standard Platform	56
SEVEN Drill Stopper Kits, Depth Based	57
Bone Compression Kit	58
Prosthetic Tools Kit	60
Abutment Holding System Kit	63



PROSTHETIC OPTIONS	64 - 93.
Analogs and Impression Copings	66
Temporary Cylinders	68
Multi-Unit System	70
Plastic Cylinders	73
EZ-Base System	75
Ti-Bases	76
Scan Posts	77
Digital Model Analogs	78
Titanium Blanks	79
Cementable Abutments	80
Cementing Posts	83
CPK - Complete Prosthetic Kit	84
LOCKiT Attachment System	87
OT-Equator Kits	90
Ball Attachment System	92



ADDITIONAL INFORMATION	106 - 116.
Pilot Drills	108
SEVEN Twist Drills	109
M4 Twist Drills	110
SEVEN and M4 Short Twist Drills	111
Alloy Composition	112



INDEX	113-120
-------	---------



MGUIDE	94- 105.
MGUIDE Templates	96
MGUIDE SEVEN Surgical Set	97
MGUIDE Drill Replacement Kit	101
MGUIDE SEVEN Surgical Set, Narrow Sleeve	102
MGUIDE Drill Replacement Kit, Narrow Sleeve	105



IMPLANT SYSTEMS

MIS offers a range of implant systems that are simple and reliable. MIS implant systems seamlessly Integrate into a digital workflow, allowing clinicians to offer their patients advanced treatment options.



The surface of all MIS implants is sand-blasted and acid etched. The surface roughness, uniformity, cleanliness of the implants are regularly monitored, with an in-house Scanning Electron Microscope.





Implant Systems

SEVEN XD Implant Internal Hex. Ø3.30 NP




	Cat. No	Dimensions (mm)	Description	Material	Item includes
	MF7-D10330	Ø3.30 L=10	SEVEN XD internal hex. implant, narrow platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills
	MF7-D11330	Ø3.30 L=11.5	SEVEN XD internal hex. implant, narrow platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills
	MF7-D13330	Ø3.30 L=13	SEVEN XD internal hex. implant, narrow platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills
	MF7-D16330	Ø3.30 L=16	SEVEN XD internal hex. implant, narrow platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills

L = Length

Ø3.30mm Drilling Protocol

Drilling Speed (RPM)	1200-1500	900-1200	200-400	15-25
Diameter	Ø1.90	Ø2.40-2.80	Ø3.25	Ø3.30





i Do not use the last drill for bone types 3 and 4. The drilling sequence is illustrated using a 13mm implant. Procedures recommended by MIS cannot replace the judgment and professional experience of the surgeon.



Implant Systems

SEVEN XD Implant Internal Hex. Ø3.75 SP




	Cat. No	Dimensions (mm)	Description	Material	Item includes
	MF7-D08375	Ø3.75 L=8	SEVEN XD internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills
	MF7-D10375	Ø3.75 L=10	SEVEN XD internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills
	MF7-D11375	Ø3.75 L=11.5	SEVEN XD internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills
	MF7-D13375	Ø3.75 L=13	SEVEN XD internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills
	MF7-D16375	Ø3.75 L=16	SEVEN XD internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills

L = Length

Ø3.75mm Drilling Protocol

Drilling Speed (RPM)	1200-1500	900-1200	200-400	200-400	15-25
Diameter	Ø1.90	Ø2.40-2.80	Ø3.25	Ø3.65	Ø3.75





Do not use the last drill for bone types 3 and 4. The drilling sequence is illustrated using a 13mm implant. Procedures recommended by MIS cannot replace the judgment and professional experience of the surgeon.



Implant Systems

SEVEN XD Implant Internal Hex. Ø4.20 SP




	Cat. No	Dimensions (mm)	Description	Material	Item includes
	MF7-D06420	Ø4.20 L=6	SEVEN XD internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills
	MF7-D08420	Ø4.20 L=8	SEVEN XD internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills
	MF7-D10420	Ø4.20 L=10	SEVEN XD internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills
	MF7-D11420	Ø4.20 L=11.5	SEVEN XD internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills
	MF7-D13420	Ø4.20 L=13	SEVEN XD internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills
	MF7-D16420	Ø4.20 L=16	SEVEN XD internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills

L = Length

Ø4.20mm Drilling Protocol

Drilling Speed (RPM)	1200-1500	900-1200	200-400	200-400	15-25
Diameter	Ø1.90	Ø2.40-2.80	Ø3.65	Ø4.10	Ø4.20




Do not use the last drill for bone types 3 and 4. The drilling sequence is illustrated using a 13mm implant. Procedures recommended by MIS cannot replace the judgment and professional experience of the surgeon.



Implant Systems

SEVEN XD Implant Internal Hex. Ø5 WP



Implant Cover Screw
and Healing Caps:

Standard



MH-W3500
MH-W4500
MH-W5500

Standard



MH-W0455
MH-W0555
MH-W0655

Anatomic



MH-W3630
MH-W4630
MH-W5630

Cover screw




MC1-00470

	Cat. No	Dimensions (mm)	Description	Material	Item includes
	MF7-D06500	Ø5 L=6	SEVEN XD internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills
	MF7-D08500	Ø5 L=8	SEVEN XD internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills
	MF7-D10500	Ø5 L=10	SEVEN XD internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills
	MF7-D11500	Ø5 L=11.5	SEVEN XD internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills
	MF7-D13500	Ø5 L=13	SEVEN XD internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills
	MF7-D16500	Ø5 L=16	SEVEN XD internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw XD single-use drills

L = Length

Ø5mm Drilling Protocol

Drilling Speed (RPM)	1200-1500	900-1200	200-400	200-400	200-400	200-400	15-25
Diameter	Ø1.90	Ø2.40-2.80	Ø3.65	Ø4.10	Ø4.9*	Ø4.9	Ø5







i * When placing Ø5 implants in soft bone, the last drill should be drilled to the first depth indicator, which is 6mm deep. For 8mm length implants, the first depth marking is 4mm deep.

Do not use the last drill for bone types 3 and 4. The drilling sequence is illustrated using a 13mm implant. Procedures recommended by MIS cannot replace the judgment and professional experience of the surgeon.



Implant Systems

SEVEN Implant Internal Hex. Ø6 WP



	Cat. No	Dimensions (mm)	Description	Material	Item includes
	MF7-06600	Ø6 L=6	SEVEN internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw Final drill
	MF7-08600	Ø6 L=8	SEVEN internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw Final drill
	MF7-10600	Ø6 L=10	SEVEN internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw Final drill
	MF7-11600	Ø6 L=11.5	SEVEN internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw Final drill
	MF7-13600	Ø6 L=13	SEVEN internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw Final drill

L = Length


* SEVEN 6mm implants are not available with XD drills.

Ø6mm Drilling Protocol

Drilling Speed (RPM)	1200- 1500	900- 1200		500- 700	400- 700	400- 600	300- 500	300- 500		200-400	200-500	15-25
Diameter	Ø1.90	Ø2.40	Ø2.40	Ø2.80	Ø3.20	Ø4	Ø4.50	Ø5	Ø5	Ø5.30 Ø5.90	Ø6	Ø6

Hard Bone (Type 1&2)

Soft Bone (Type 3&4)



For SEVEN Ø6mm implants, use the SEVEN surgical kit (MK-T048).

Do not use the final drill for bone types 3 and 4. The drilling sequence is illustrated using a 13mm implant. Procedures recommended by MIS cannot replace the judgment and professional experience of the surgeon.



Implant Systems

XD Single-Use Drills

REVOLUTION IN A BOX!

Unique in the market, MIS provides sharp and sterile XD drills with each SEVEN implant, designed for single use.

SINGLE-USE

XD drills are used once per procedure, ensuring strict clinical hygiene protocols for clinicians and patients. Using single-use drills may reduce the risk of possible cross-contamination and infection.

SHARP

Using sharp drills in every drilling procedure prevents drill wear and deformation.

STERILE

XD drills eliminate the need for pre-surgery sterilization.

SOLID

XD drills are designed for depth control, with clear depth markings, providing visibility and promoting confidence in every drilling procedure.

SIMPLE

The MIS XD drilling procedure is simplified and quick, eliminating the need for separate cleaning and sterilization in a surgical kit, as the drills are sterile and for single use.



MIS XD GUARANTEED QUALITY

MIS has performed bench testing of XD drills, including:

- Cutting ability test: MIS XD have shown cutting ability similar to reusable kit drills, according to internal testing procedures.
- Corrosion test: MIS XD have shown no corrosion throughout its shelf-life.
- Bending test: MIS XD meet the requirements of ISO 8325:2004.
- Drilling ability test: 30 clinicians evaluated the drills in simulated use; 29 of them have reported favorable results.
- Hardness test: MIS XD meet the requirements of ISO 1797:2017.

Bone Protection

Studies [1-3] have shown that bone temperature increase is directly correlated with higher drill wear. MIS XD are new and sharp each time.

Depth control

XD procedure drills are 16mm long and marked to indicate the MIS implant lengths: ,11.5 ,10 ,8 ,6 13 and 16mm.

XD final drill

Correspond to the implant's shape and dimensions.

Precise Positioning





The tip of the XD drill has a step tip with the same diameter as the previously-used drill.

[1] Chacon GE, Bower DL, Larsen PE, McGlumphy EA, Beck FM. Heat production by 3 implant drill systems after repeated drilling and sterilization. J Oral Maxillofac Surg. 2006 Feb;64(2):265-269.
[2] Misir AF, Sumer M, Yenisey M, Ergioglu E. Effect of surgical drill guide on heat generated from implant drilling. J Oral Maxillofac Surg. 2009 Dec;67(12):2663-2668.
[3] Möhlhenrich SC, Modabber A, Steiner T, Mitchell DA, Hölzle F. Heat generation and drill wear during dental implant site preparation: systematic review. Br J Oral Maxillofac Surg. 2015 Oct;53(8):679-689.



Implant Systems

M4 Implant Internal Hex. Ø3.30 NP




	Cat. No	Dimensions (mm)	Description	Material	Item includes
	MF4-10330	Ø3.30 L=10	M4 internal hex. implant, narrow platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-11330	Ø3.30 L=11.5	M4 internal hex. implant, narrow platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-13330	Ø3.30 L=13	M4 internal hex. implant, narrow platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-16330	Ø3.30 L=16	M4 internal hex. implant, narrow platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw

L = Length

Ø3.30mm Drilling Protocol



Drilling Speed (RPM)	1200- 1500	900- 1200	500-700	200-500	15-25
Diameter	Ø1.90	Ø2.40	Ø2.80	Ø3 Ø3.30	Ø3.30





Implant Systems

M4 Implant
Internal Hex.
Ø3.75 SP



Implant Cover Screw and Healing Caps:


Thin	Standard	Standard	Anatomic	Cover screw
 Ø4mm 3/4/5/6mm	 Ø4.8mm 3/4/5/6mm	 Ø4.8mm 4/5/6mm CONCAVE	 Ø5.5mm 3/4/5/6mm	
MH-03375 MH-04375 MH-05375 MH-06375	MH-A3375 MH-A4375 MH-A5375 MH-A6375	MH-S0448 MH-S0548 MH-S0648	MH-53375 MH-54375 MH-55375 MH-56375	MC1-00375

	Cat. No	Dimensions (mm)	Description	Material	Item includes
	MF4-08375	Ø3.75 L=8	M4 internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-10375	Ø3.75 L=10	M4 internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-11375	Ø3.75 L=11.5	M4 internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-13375	Ø3.75 L=13	M4 internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-16375	Ø3.75 L=16	M4 internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw

L = Length

Ø3.75mm Drilling Protocol

Drilling Speed (RPM)	1200- 1500	900- 1200	500- 700	400-700	200-500	15-25
Diameter	Ø1.90	Ø2.40	Ø2.80	Ø3.20	Ø3.75	Ø3.75





i The drilling sequence is illustrated using a 13mm implant.
Procedure recommended by MIS cannot replace the judgment and professional experience of the surgeon.



Implant Systems

M4 Implant
Internal Hex.
Ø4.20 SP



Implant Cover Screw and Healing Caps:


Thin	Standard	Standard	Anatomic	Cover screw
 Ø4mm 3/4/5/6mm	 Ø4.8mm 3/4/5/6mm	 Ø4.8mm 4/5/6mm CONCAVE	 Ø5.5mm 3/4/5/6mm	
MH-03375 MH-04375 MH-05375 MH-06375	MH-A3375 MH-A4375 MH-A5375 MH-A6375	MH-S0448 MH-S0548 MH-S0648	MH-53375 MH-54375 MH-55375 MH-56375	MC1-00375

	Cat. No	Dimensions (mm)	Description	Material	Item includes
	MF4-06420	Ø4.20 L=6	M4 internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-08420	Ø4.20 L=8	M4 internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-10420	Ø4.20 L=10	M4 internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-11420	Ø4.20 L=11.5	M4 internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-13420	Ø4.20 L=13	M4 internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-16420	Ø4.20 L=16	M4 internal hex. implant, standard platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw

L = Length

Ø4.20mm Drilling Protocol

Drilling Speed (RPM)	1200- 1500	900- 1200	500- 700	400- 700	400-600	200-500	15-25
Diameter	Ø1.90	Ø2.40	Ø2.80	Ø3.20	Ø3.80	Ø4.20	Ø4.20







 The drilling sequence is illustrated using a 13mm implant.
Procedure recommended by MIS cannot replace the judgment and professional experience of the surgeon.



Implant Systems

M4 Implant Internal Hex. Ø5 WP




	Cat. No	Dimensions (mm)	Description	Material	Item includes
	MF4-06500	Ø5 L=6	M4 internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-08500	Ø5 L=8	M4 internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-10500	Ø5 L=10	M4 internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-11500	Ø5 L=11.5	M4 internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-13500	Ø5 L=13	M4 internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-16500	Ø5 L=16	M4 internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw

L = Length

Ø5mm Drilling Protocol

Drilling Speed (RPM)	1200-1500	900-1200	500-700	400-700	400-600	400-600	200-500	15-25
Diameter	Ø1.90	Ø2.40	Ø2.80	Ø3.20	Ø3.80	Ø4.50	Ø5	Ø5





i The drilling sequence is illustrated using a 13mm implant.
Procedure recommended by MIS cannot replace the judgment and professional experience of the surgeon.



Implant Systems

M4 Implant
Internal Hex.
Ø6 WP




	Standard	Standard	Anatomic	Cover screw
Implant Cover Screw and Healing Caps:				
	MH-W3500 MH-W4500 MH-W5500	MH-W0455 MH-W0555 MH-W0655	MH-W3630 MH-W4630 MH-W5630	MC1-00470

	Cat. No	Dimensions (mm)	Description	Material	Item includes
	MF4-06600	Ø6 L=6	M4 internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-08600	Ø6 L=8	M4 internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-10600	Ø6 L=10	M4 internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-11600	Ø6 L=11.5	M4 internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw
	MF4-13600	Ø6 L=13	M4 internal hex. implant, wide platform	Titanium Alloy Ti 6Al 4V ELI	Cover screw

L = Length

Ø6mm Drilling Protocol

Drilling Speed (RPM)	1200-1500	900-1200	500-700	400-700	400-600	400-600	400-600	300-500	200-500	15-25
Diameter	Ø1.90	Ø2.40	Ø2.80	Ø3.20	Ø3.80	Ø4.50	Ø5	Ø5.50	Ø6	Ø6



The drilling sequence is illustrated using a 13mm implant.
Procedure recommended by MIS cannot replace the judgment and professional experience of the surgeon.



TOOLS AND KITS

MIS kits and instruments promote accuracy and precision in surgical procedures. The kits are designed for efficiency and ease of use, featuring an ergonomic design and a single-handed opening mechanism. The range includes various tools such as drills, body try-ins, keys, ratchets, and insertion tools.





Tools and Kits

Step Pilot Drills




	Cat. No	Dimensions (mm)	Description	Material	Sold
	CT-P2406	Ø2.4/2.0 L=31.8	Step pilot drill with built-in stopper for 6mm length implants	Stainless steel	 
	CT-P2408	Ø2.4/2.0 L=28	Step pilot drill with built-in stopper for 8mm length implants	Stainless steel	 
	CT-P2410	Ø2.4/2.0 L=28	Step pilot drill with built-in stopper for 10mm length implants	Stainless steel	 

L = Length

	Cat. No	Dimensions (mm)	Description	Material	Sold
	CT-P2411	Ø2.4/2.0 L=33	Step pilot drill with built-in stopper for 11.5mm length implants	Stainless steel	 
	CT-P2413	Ø2.4/2.0 L=33	Step pilot drill with built-in stopper for 13mm length implants	Stainless steel	 
	CT-P2416	Ø2.4/2.0 L=37.5	Step pilot drill for 16mm length implants	Stainless steel	 

L = Length

6 Depth Marks





Tools and Kits

SEVEN Twist Drills



Cat. No	Dimensions (mm)	Description	Material	Sold
MT-TDT28	Ø2.80 L=37.5	Twist drill	Stainless steel	 
MT-TDT32	Ø3.20 L=37.5	Twist drill	Stainless steel	 
MT-TDT40	Ø4 L=37.5	Twist drill	Stainless steel	 
MT-TDT45	Ø4.50 L=38.2	Twist drill	Stainless steel	 
MT-TDT50	Ø5 L=38.2	Twist drill	Stainless steel	 

L = Length



Tools and Kits

M4 Twist Drills



Cat. No	Dimensions (mm)	Description	Material	Sold
MT-TDN28	Ø2.80 L=37.5	Twist drill	Stainless steel	 
MT-TDN32	Ø3.20 L=37.6	Twist drill	Stainless steel	 
MT-TDN38	Ø3.80 L=37.8	Twist drill	Stainless steel	 
MT-TDN45	Ø4.50 L=38.2	Twist drill	Stainless steel	 
MT-TDN50	Ø5 L=38.2	Twist drill	Stainless steel	 
MT-TDN55	Ø5.50 L=38.2	Twist drill	Stainless steel	 

L = Length



Tools and Kits

SEVEN and M4 Short Twist Drills




	Cat. No	Dimensions (mm)	Description	Material	Sold
	MT-TDI28	Ø2.80 L=30	Short twist drill	Stainless steel	
	MT-TDI32	Ø3.20 L=30	Short twist drill	Stainless steel	
	MT-TDI38	Ø3.80 L=30	Short twist drill	Stainless steel	
	MT-TDI40	Ø4.0 L=30	Short twist drill	Stainless steel	

L = Length

	Cat. No	Dimensions (mm)	Description	Material	Sold
	MT-TDI45	Ø4.50 L=30	Short twist drill	Stainless steel	
	MT-TDI50	Ø5.0 L=30	Short twist drill	Stainless steel	
	MT-TDI55	Ø5.50 L=30	Short twist drill	Stainless steel	

L = Length

7 Drilling Diameter Options





Tools and Kits

Marking Drills and Countersinks



	Cat. No	Dimensions (mm)	Description	Material	Sold
	MT-TDN19	Ø1.90 L=34	Marking drill	Stainless steel	
	MT-SMD10	L=30	Spade marking drill	Stainless steel	
	MT-PD440	Ø4 L=32.7	Position drill	Stainless steel	

L = Length


Cat. No	Dimensions (mm)	Description	Material	Sold
MT-CS330	Ø3.30 L=25	Countersink for Ø3.30 implant, narrow platform	Stainless steel	STAND ALONE
MT-CS375	Ø3.75 L=25	Countersink for Ø3.75 implant, standard platform	Stainless steel	STAND ALONE
MT-CS420	Ø4.20 L=25	Countersink for Ø4.20 implant, standard platform	Stainless steel	STAND ALONE
MT-CS500	Ø5 L=25	Countersink for Ø5 implant, wide platform	Stainless steel	STAND ALONE
MT-CSN33	Ø3.30 L=31	Countersink, narrow platform	Stainless steel	INSIDE KIT STAND ALONE
MT-GDN33	Ø3.75/4.20 L=26	Countersink, standard platform	Stainless steel	INSIDE KIT STAND ALONE
MT-GDN50	Ø5/6 L=26.8	Countersink, wide platform	Stainless steel	INSIDE KIT STAND ALONE

L = Length




5



Countersinks Depth Marks



Narrow Platform



Standard Platform



Wide Platform



Tools and Kits

Trephine Burrs



Cat. No	Dimensions (mm)	Description	Material	Sold
---------	-----------------	-------------	----------	------

MT-00400	Ø4/5 L=37	Trephine burr	Stainless steel	STAND ALONE
MT-00500	Ø5/6 L=37	Trephine burr	Stainless steel	STAND ALONE
MT-00600	Ø6/7 L=37	Trephine burr	Stainless steel	STAND ALONE



L = Length



Tools and Kits

Insertion Tools



	Cat. No	Dimensions (mm)	Description	Material	Sold
	MT-HLI21	L=28.2	Long motor insertion tool, internal hex., narrow platform	Stainless steel	 
	MT-HSI21	L=24.4	Short motor insertion tool, internal hex., narrow platform	Stainless steel	 
	MT-HLI10	L=28.2	Long motor insertion tool, internal hex., standard/wide platform	Stainless steel	 
	MT-HSI10	L=24.4	Short motor insertion tool, internal hex., standard/wide platform	Stainless steel	 
	MT-LRH21	L=32.3	Long ratchet insertion tool, internal hex., narrow platform	Stainless steel	 
	MT-SRH21	L=22.3	Short ratchet insertion tool, internal hex., narrow platform	Stainless steel	 
	MT-LRH20	L=32.3	Long ratchet insertion tool, internal hex., standard/wide platform	Stainless steel	 
	MT-SRH20	L=22.3	Short ratchet insertion tool, internal hex., standard/wide platform	Stainless steel	 

L = Length



Tools and Kits

Surgical Tools



Cat. No	Dimensions (mm)	Description	Material	Sold
MT-BTI20	Ø1.40 L=100	Implant site depth probe	Stainless steel	INSIDE KIT STAND ALONE
MT-RDE30	L=32	Extra-long driver for 0.05 inch hex.	Stainless steel	STAND ALONE
MT-RDL30	L=27	Long driver for 0.05 inch hex.	Stainless steel	INSIDE KIT STAND ALONE
MT-RDS30	L=21.3	Short driver for 0.05 inch hex.	Stainless steel	INSIDE KIT STAND ALONE

L = Length

	Cat. No	Dimensions (mm)	Description	Material	Sold
           	MT-RDM30	L=11.5	Extra-short driver for 0.05 inch hex.	Stainless steel	
	MT-DE001	L=24	Drill extender	Stainless steel	 
	MT-PP240	Ø2.40-3 L=24	Parallel pin	Stainless steel	
	MT-LM005	L=30	Long motor adapter for 0.05 inch hex.	Stainless steel	 
	MT-SM005	L=24	Short motor adapter for 0.05 inch hex.	Stainless steel	
	MT-SRA10	L=15.5	Square connection to ratchet adapter	Stainless steel	
	MT-TP420	Ø4.30 L=30	Motor tissue punch	Stainless steel	
	MT-TP500	Ø4.30 L=28	Manual tissue punch	Stainless steel	
	MT-RT070	L=95	Surgical torque ratchet		 
	MT-RI030	L=84	Monoblock ratchet		 
	MT-RA480	L=7.5	Ratchet to hex. adapter	Stainless steel	
	MT-MFXGS		MULTIFIX guide strip	Stainless steel	

L = Length



Tools and Kits

Prosthetic Instruments



	Cat. No	Dimensions (mm)	Description	Material	Sold
	MT-IT100	L=16.5	Impression coping key for closed tray	Stainless steel	
	MT-LM005	L=30	Long motor adapter for 0.05 inch hex.	Stainless steel	 
	MT-SM005	L=24	Short motor adapter for 0.05 inch hex.	Stainless steel	
	MT-MSD20	L=140	Surgical adapter driver	Stainless steel + Aluminum	 

L = Length

	Cat. No	Dimensions (mm)	Description	Material	Sold
	MT-UN002	L=24	Implant analog for universal hand driver, internal hex., narrow platform	Stainless steel	 
	MT-US002	L=24	Implant analog for universal hand driver, internal hex., standard platform	Stainless steel	 
	MT-UW002	L=24	Implant analog for universal hand driver, internal hex., wide platform	Stainless steel	 
	MT-HTS10	L=144,5	Technical hex. screwdriver	Stainless steel + plastic	 
	MT-RB225	L=22	Ball attachment anchor key	Stainless steel	 
	MT-RDE30	L=32	Extra-long driver for 0.05 inch hex.	Stainless steel	
	MT-RDL30	L=27	Long driver for 0.05 inch hex.	Stainless steel	 
	MT-RDS30	L=21,3	Short driver for 0.05 inch hex.	Stainless steel	 
	MT-RDM30	L=11,5	Extra-short driver for 0.05 inch hex.	Stainless steel	
	MT-MURL2	L=21	Straight Multi-Unit long ratchet key	Stainless steel	 
	MT-MURS2	L=16	Straight Multi-Unit short ratchet key	Stainless steel	
	MT-MUML2	L=31	Straight Multi-Unit long motor key	Stainless steel	

L = Length

Prosthetic Instruments


	Cat. No	Dimensions (mm)	Description	Material	Sold
	MT-MUMS2	L=24	Straight Multi-Unit short motor key	Stainless steel	STAND ALONE
	MT-RSD19	L=19	LOCKiT short torque ratchet driver	Stainless steel	INSIDE KIT STAND ALONE
	MT-RSD24	L=24	LOCKiT long torque ratchet driver	Stainless steel	INSIDE KIT STAND ALONE
	MT-ELM10	L=30	Long motor key for angled screw channel	Stainless steel	STAND ALONE
	MT-ESM10	L=24	Short motor key for angled screw channel	Stainless steel	STAND ALONE
	MT-ELR10	L=22.5	Long ratchet key for angled screw channel	Stainless steel	INSIDE KIT STAND ALONE
	MT-ESR10	L=18.7	Short ratchet key for angled screw channel	Stainless steel	STAND ALONE
	MT-LOC10	Ø7.8 L=86	LOCKiT core tool	Stainless steel	INSIDE KIT STAND ALONE
	ET-IT002		OT-Equator handling tool		INSIDE KIT STAND ALONE
	MG1-TE480	Ø8.6 L=16	MGUIDE tool extender	Stainless steel	STAND ALONE
	MT-RA480		Ratchet to hex. adapter	Stainless steel	STAND ALONE
	MT-RI040	L=95	Torque ratchet for prosthetic screws		INSIDE KIT STAND ALONE

L = Length



Tools and Kits

SEVEN XD Placement Set



MK-XD01
XD organizer tray

The XD Organizer Tray is designed to support the XD procedure by providing a convenient place to store the XD drills during surgery. Please note that the organizer isn't included in the surgical kit and should be sterilized prior to surgery.

MK-SVXD
SEVEN XD Placement Set



Cat. No.	Dimensions (mm)	Description	Material	Sold
MT-PDM24	Ø2.4 L=32	Position drill mill	Stainless steel	 
MT-TDN19	Ø1.90 L=34	Marking drill	Stainless steel	 
CT-P2416	Ø2.4/2.0 L=37.5	Step pilot drill for 16mm length implants	Stainless steel	 
CT-BTC24	Ø2.4 L=28.5	Body try-in	Titanium	 

L = Length

SEVEN XD Placement Set

	Cat. No.	Dimensions (mm)	Description	Material	Sold
	MT-CSN33	Ø3.3 L=31	Countersink, narrow platform	Stainless steel	 
	MT-GDN33	Ø3.75/4.2 L=26	Countersink, standard platform	Stainless steel	 
	MT-GDN50	Ø5/6 L=26.8	Countersink, wide platform	Stainless steel	 
	MT-DE001	L=24	Drill extender	Stainless steel	 
	MT-RDL30	L=27	Long driver for 0.05 inch hex.	Stainless steel	 
	MT-RDS30	L=21.3	Short driver for 0.05 inch hex.	Stainless steel	 
	MT-HLI21	L=28.2	Long motor insertion tool, internal hex., narrow platform	Stainless steel	 
	MT-HLI10	L=28.2	Long motor insertion tool, internal hex., standard/wide platform	Stainless steel	 
	MT-LRH21	L=32.3	Long ratchet insertion tool, internal hex., narrow platform	Stainless steel	 
	MT-LRH20	L=32.3	Long ratchet insertion tool, internal hex., standard/wide platform	Stainless steel	 
	MT-RMR10	L=38.5	Long direct hand and ratchet key, standard/wide platform	Stainless steel	 
	MT-RT070	L=95	Surgical torque ratchet		 

L = Length



Tools and Kits

SEVEN Surgical Kit

MK-T048
SEVEN Surgical Kit



Cat. No	Dimensions (mm)	Description	Material	Sold
MT-TDN19	Ø1.90 L=34	Marking drill	Stainless steel	INSIDE KIT STAND ALONE
MT-SMD10	L=30	Spade marking drill	Stainless steel	INSIDE KIT STAND ALONE
CT-P2406	Ø2.40/2.0 L=31.8	Step pilot drill with built-in stopper for 6mm length implants	Stainless steel	INSIDE KIT STAND ALONE
CT-P2408	Ø2.40/2.0 L=28	Step pilot drill with built-in stopper for 8mm length implants	Stainless steel	INSIDE KIT STAND ALONE
CT-P2410	Ø2.40/2.0 L=28	Step pilot drill with built-in stopper for 10mm length implants	Stainless steel	INSIDE KIT STAND ALONE
CT-P2411	Ø2.40/2.0 L=33	Step pilot drill with built-in stopper for 11.5mm length implants	Stainless steel	INSIDE KIT STAND ALONE

L = Length

SEVEN
Surgical Kit

	Cat. No	Dimensions (mm)	Description	Material	Sold
	CT-P2413	Ø2.40/2.0 L=33	Step pilot drill with built-in stopper for 13mm length implants	Stainless steel	 
	CT-P2416	Ø2.40/2.0 L=37.5	Step pilot drill for 16mm length implants	Stainless steel	 
	MT-TDT28	Ø2.80 L=37.5	Twist drill	Stainless steel	 
	MT-TDT32	Ø3.20 L=37.5	Twist drill	Stainless steel	 
	MT-TDT40	Ø4 L=37.5	Twist drill	Stainless steel	 
	MT-TDT45	Ø4.50 L=38.2	Twist drill	Stainless steel	 
	MT-TDT50	Ø5 L=38.2	Twist drill	Stainless steel	 
	CT-BTC24	Ø2.40 L=28.5	Body try-in	Titanium	 
	MT-BTT28	Ø2.80 L=28.5	Body try-in	Titanium	 
	MT-BTT32	Ø3.20 L=28.5	Body try-in	Titanium	 
	MT-BTT40	Ø4 L=28.5	Body try-in	Titanium	 
	MT-BTT45	Ø4.50 L=28.5	Body try-in	Titanium	 

L = Length

	Cat. No	Dimensions (mm)	Description	Material	Sold
	MT-BTT50	Ø5 L=28.5	Body try-in	Titanium	 
	MT-GDN33	Ø3.75/4.20 L=26	Countersink, standard platform	Stainless steel	 
	MT-GDN50	Ø5/6 L=26.8	Countersink, wide platform	Stainless steel	 
	MT-PP240	Ø2.40-3 L=24	Parallel pin	Stainless steel	 
	MT-HLI21	L=28.2	Long motor insertion tool, internal hex., narrow platform	Stainless steel	 
	MT-HLI10	L=28.2	Long motor insertion tool, internal hex., standard/wide platform	Stainless steel	 
	MT-HSI10	L=24.4	Short motor insertion tool, internal hex., standard/wide platform	Stainless steel	 
	MT-LRH21	L=32.3	Long ratchet insertion tool, internal hex., narrow platform	Stainless steel	 
	MT-LRH20	L=32.3	Long ratchet insertion tool, internal hex., standard/wide platform	Stainless steel	 
	MT-SRH20	L=22.3	Short ratchet insertion tool, internal hex., standard/wide platform	Stainless steel	 
	MT-DE001	L=24	Drill extender	Stainless steel	 
	MT-RMR10	L=38.5	Long direct hand and ratchet key, standard/wide platform	Stainless steel	 


L = Length

SEVEN
Surgical Kit

	Cat. No	Dimensions (mm)	Description	Material	Sold
	MT-RDL30	L=27	Long driver for 0.05 inch hex.	Stainless steel	 
	MT-RDS30	L=21.3	Short driver for 0.05 inch hex.	Stainless steel	 
	MN-PF330 MD-PF375	Ø4.5, L=11 Ø5, L=11	Press fit for closed tray, internal hex.	Titanium	 
	MT-RT070	L=95	Surgical torque ratchet		 

L = Length

 5 Drilling Procedures





Tools and Kits

SEVEN Drill Replacement Kit



MK-SD01
SEVEN Drill Replacement Kit

	Cat. No	Dimensions (mm)	Description	Material	Sold
	CT-P2408	Ø2.40/2.0 L=28	Step pilot drill with built-in stopper for 8mm length implants	Stainless steel	 
	CT-P2410	Ø2.40/2.0 L=28	Step pilot drill with built-in stopper for 10mm length implants	Stainless steel	 
	CT-P2411	Ø2.40/2.0 L=33	Step pilot drill with built-in stopper for 11.5mm length implants	Stainless steel	 
	CT-P2413	Ø2.40/2.0 L=33	Step pilot drill with built-in stopper for 13mm length implants	Stainless steel	 
	CT-P2416	Ø2.40/2.0 L=37.5	Step pilot drill for 16mm length implants	Stainless steel	 
	MT-TDT28	Ø2.80 L=37.5	Twist drill	Stainless steel	 
	MT-TDT32	Ø3.20 L=37.5	Twist drill	Stainless steel	 
	MT-TDT40	Ø4 L=37.5	Twist drill	Stainless steel	 
	MT-TDT45	Ø4.50 L=38.2	Twist drill	Stainless steel	 

L = Length



Tools and Kits

M4 Surgical Kit



Cat. No	Dimensions (mm)	Description	Material	Sold
MT-TDN19	Ø1.90 L=34	Marking drill	Stainless steel	INSIDE KIT STAND ALONE
MT-SMD10	L=30	Spade marking drill	Stainless steel	INSIDE KIT STAND ALONE
CT-P2416	Ø2.40/2.0 L=37.5	Step pilot drill for 16mm length implants	Stainless steel	INSIDE KIT STAND ALONE
MT-TDN28	Ø2.80 L=37.5	Twist drill	Stainless steel	INSIDE KIT STAND ALONE
MT-TDN32	Ø3.20 L=37.6	Twist drill	Stainless steel	INSIDE KIT STAND ALONE

L = Length



Cat. No	Dimensions (mm)	Description	Material	Sold
MT-TDN38	Ø3.80 L=37.8	Twist drill	Stainless steel	INSIDE KIT STAND ALONE
MT-TDN45	Ø4.50 L=38.2	Twist drill	Stainless steel	INSIDE KIT STAND ALONE
MT-TDN50	Ø5 L=38.2	Twist drill	Stainless steel	INSIDE KIT STAND ALONE
MT-TDN55	Ø5.50 L=38.2	Twist drill	Stainless steel	INSIDE KIT STAND ALONE
MT-DE001	L=24	Drill extender	Stainless steel	INSIDE KIT STAND ALONE
MT-LM005	L=30	Long motor adapter for 0.05 inch hex.	Stainless steel	INSIDE KIT STAND ALONE
MT-HSI10	L=24.4	Short motor insertion tool, internal hex., standard/wide platform	Stainless steel	INSIDE KIT STAND ALONE
MT-HLI10	L=28.2	Long motor insertion tool, internal hex., standard/wide platform	Stainless steel	INSIDE KIT STAND ALONE
MT-CSN33	Ø3.30 L=31	Countersink, narrow platform	Stainless steel	INSIDE KIT STAND ALONE
MT-GDN33	Ø3.75/4.20 L=26	Countersink, standard platform	Stainless steel	INSIDE KIT STAND ALONE
MT-GDN50	Ø5/6 L=26.8	Countersink, wide platform	Stainless steel	INSIDE KIT STAND ALONE
MT-PP240	Ø2.40-3 L=24	Parallel pin	Stainless steel	STAND ALONE
MT-LRH21	L=32.3	Long ratchet insertion tool, internal hex., narrow platform	Stainless steel	INSIDE KIT STAND ALONE

L = Length


M4 Surgical Kit

Cat. No	Dimensions (mm)	Description	Material	Sold
 MT-LRH20	L=32.3	Long ratchet insertion tool, standard/wide platform, internal hex.	Stainless steel	 
 MT-SRH20	L=22.3	Short ratchet insertion tool, standard/wide platform, internal hex.	Stainless steel	 
 MT-RMR10	L=38.5	Long direct hand and ratchet key, standard/wide platform	Stainless steel	 
 MT-RDL30	L=27	Long driver for 0.05 inch hex.	Stainless steel	 
 MT-RI030	L=84	Monoblock ratchet		 

L = Length



5 Surgical Procedures





Tools and Kits

Bone Profiler Kit

MK-0059
Bone Profiler Kit



	Cat. No	Dimensions (mm)	Description	Material	Sold
	MT-BP320	Ø3.20	Bone profiler, internal hex., narrow platform	Stainless steel	INSIDE KIT
	MT-BP370	Ø3.70	Bone profiler, internal hex., standard platform	Stainless steel	INSIDE KIT
	MT-BP460	Ø4.60	Bone profiler, internal hex., wide platform	Stainless steel	INSIDE KIT
	MT-BGP32	Ø3.20	Bone profiler guide pin, internal hex., narrow platform	Stainless steel	INSIDE KIT
	MT-BGP37	Ø3.70	Bone profiler guide pin, internal hex., standard platform	Stainless steel	INSIDE KIT
	MT-BGP46	Ø4.60	Bone profiler guide pin, internal hex., wide platform	Stainless steel	INSIDE KIT

L = Length




Tools and Kits

SEVEN Drill Stoppers Kit, Standard Platform

Following the use of MIS length-based pilot drills, MIS Drill Stopper Kits ensure drilling to the desired depth.

The Standard Platform Kit includes the most commonly used implant lengths: 8, 10, 11.5 and 13mm.



Cat. No

Description



MK-BS001

Drill stopper kit for SEVEN implant, standard platform




Tools and Kits

SEVEN Drill Stopper Kits, Depth Based

MK-SDS(06-13) 
Drill Stopper Kits



	Cat. No	Dimensions (mm)	Description	Material
	MK-SDS06	L=6	Drill stoppers kit for 6mm SEVEN implant	Stainless steel
	MK-SDS08	L=8	Drill stoppers kit for 8mm SEVEN implant	Stainless steel
	MK-SDS10	L=10	Drill stoppers kit for 10mm SEVEN implant	Stainless steel
	MK-SDS11	L=11.5	Drill stoppers kit for 11.5mm SEVEN implant	Stainless steel
	MK-SDS13	L=13	Drill stoppers kit for 13mm SEVEN implant	Stainless steel

L = Length



Tools and Kits

Bone Compression Kit

MK-0020
Bone Compression Kit




	Cat. No	Dimensions (mm)	Description	Material	Sold
	MO-CO160	Ø1.80-2.90 L=25	Convex compression screw	Stainless steel	
	MO-CO200	Ø2.30-3.30 L=25	Convex compression screw	Stainless steel	
	MO-CO240	Ø2.70-3.70 L=25	Convex compression screw	Stainless steel	
	MO-CO280	Ø3.10-4.50 L=25	Convex compression screw	Stainless steel	
	MO-SO300	Ø3 L=25	Concave sinus screw	Stainless steel	
	MO-SO350	Ø3.50 L=25	Concave sinus screw	Stainless steel	
	MO-SO400	Ø4 L=25	Concave sinus screw	Stainless steel	

L = Length


Cat. No	Dimensions (mm)	Description	Material	Sold
MT-SMD10	L=30	Spade marking drill	Stainless steel	INSIDE KIT STAND ALONE
MT-HSI10	L=24.4	Short motor insertion tool, internal hex., standard/wide platform	Stainless steel	INSIDE KIT STAND ALONE
MT-HLI10	L=28.2	Long motor insertion tool, internal hex., standard/wide platform	Stainless steel	INSIDE KIT STAND ALONE
MT-LRH20	L=32.3	Long ratchet insertion tool, internal hex., standard/wide platform	Stainless steel	INSIDE KIT STAND ALONE
MT-SRH20	L=22.3	Short ratchet insertion tool, internal hex., standard/wide platform	Stainless steel	INSIDE KIT STAND ALONE
MT-RI030	L=84	Monoblock ratchet		INSIDE KIT STAND ALONE

L = Length

2 Surgical Procedures



Bone Compression



Sinus Lifting



Tools and Kits

Prosthetic Tools Kit

MK-0039
Prosthetic Tools Kit



Cat. No.	Dimensions (mm)	Description	Material	Sold
MT-RDS30	L=21.3	Short driver for 0.05 inch hex.	Stainless steel	INSIDE KIT STAND ALONE
MT-RDL30	L=27	Long driver for 0.05 inch hex.	Stainless steel	INSIDE KIT STAND ALONE
MT-SM005	L=24	Short motor adapter for 0.05 inch hex.	Stainless steel	INSIDE KIT STAND ALONE
MT-LM005	L=30	Long motor adapter for 0.05 inch hex.	Stainless steel	INSIDE KIT STAND ALONE
MT-CLR21	L=31.6	CONNECT long ratchet insertion tool	Stainless steel	INSIDE KIT STAND ALONE

L = Length



Cat. No.	Dimensions (mm)	Description	Material	Sold
MT-MURL2	L=21	Straight Multi-Unit long ratchet key	Stainless steel	 
MT-MUML2	L=31	Straight Multi-Unit long motor key	Stainless steel	 
MT-RE160	L=28.5	Abutment extractor, conical connection, narrow platform	Titanium	 
MT-RE172	L=28.5	Abutment extractor, conical connection, standard/wide platform	Titanium	 
MT-RSD19	L=19	LOCKiT short torque ratchet driver	Stainless steel	 
MT-RSD24	L=24	LOCKiT long torque ratchet driver	Stainless steel	 
MT-RB225	L=22	Ball attachment anchor key	Stainless steel	 
MT-UN002	L=24	Implant analog for universal hand driver, internal hex., narrow platform	Stainless steel	 
MT-US002	L=24	Implant analog for universal hand driver, internal hex., standard platform	Stainless steel	 
MT-UW002	L=24	Implant analog for universal hand driver, internal hex., wide platform	Stainless steel	 
CT-UN001	L=24	Implant analog for universal hand driver, conical connection, narrow platform	Stainless steel	
CT-US001	L=24	Implant analog for universal hand driver, conical connection, standard platform	Stainless steel	

L = Length

Prosthetic
Tools Kit


	Cat. No.	Dimensions (mm)	Description	Material	Sold
	CT-UW001	L=24	Implant analog for universal hand driver, conical connection, wide platform	Stainless steel	INSIDE KIT
	ET-IT002		OT-Equator handling tool		INSIDE KIT STAND ALONE
	MT-LOC10	Ø7.8 L=86	LOCKiT core tool	Stainless steel	INSIDE KIT STAND ALONE
	MT-MSD20	L=140	Surgical adapter driver	Stainless steel + Aluminum	INSIDE KIT STAND ALONE
	MT-ELR10	L=22.5	Long ratchet key for angled screw channel	Stainless steel	INSIDE KIT STAND ALONE
	MT-RI040	L=95	Torque ratchet for prosthetic screws		INSIDE KIT STAND ALONE

L = Length



4 Prosthetic Tools

2 Implant Connections





Tools and Kits

Abutment Holding System Kit

The abutment holding system allows convenient grip of prosthetic components for adjustments and preparation of both ceramic and titanium abutments.



Cat. No.

Description



MK-0001

Abutment holding system kit, internal hex.



PROSTHETIC OPTIONS

MIS prosthetic components are suitable for a wide range of clinical applications, including implant level or abutment level, one-stage or two-stage treatments, and single or multiple-unit solutions. They come in various gingival heights, diameters, and angles.

The MIS prosthetic line features a consistent concave emergence profile, designed to give more space for soft tissue ingrowth. The gold-shaded abutments minimize the reflection through the gingiva^[4-5].

[4] Ioannidis A, Cathomen E, Jung RE, Fehmer V, Hüsler J, Thoma DS. Discoloration of the mucosa caused by different restorative materials - a spectrophotometric in vitro study. Clin Oral Implants Res. 2017 Sep;28(9):1133-1138.

[5] Pecnik CM, Roos M, Muff D, Spolenak R, Sailer I. In vitro color evaluation of esthetic coatings for metallic dental implants and implant prosthetic appliances. Clin Oral Implants Res. 2015 May;26(5):563-571.





Impression Coping

Analogs and Impression Copings



	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MN-ITC30	NP	Ø4, H=2, C=6	Impression coping for closed tray, concave, internal hex.	Titanium	  
	MD-ITC30	SP	Ø4.8, H=2, C=6			
	MW-ITC30	WP	Ø5.5, H=2, C=6			
	MN-IO401	NP	Ø3.2, L=11	Impression coping for open tray, concave, internal hex.	Titanium	  
	MD-IO481	SP	Ø3.75, L=11			
	MW-IO551	WP	Ø5, L=11			
	MN-IT100	NP	Ø4, L=9	Impression coping for closed tray, internal hex.	Stainless steel	  
	MD-IT100	SP	Ø4.5, L=9			
	MW-IT100	WP	Ø5.4, L=9			
	MN-IO330	NP	Ø4.6, L=11	Impression coping for open tray, internal hex.	Stainless steel	  
	MD-IO375	SP	Ø5, L=11			
	MW-IO470	WP	Ø5.5, L=11			

 Concave emergence profile

H = Gingival height

C = Crown height

L = Length

	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MN-IT300	NP	Ø4, H=2, C=6	Impression coping set for closed tray, internal hex.	Titanium	  
	MD-IT300	SP	Ø4.8, H=2, C=6			
	MW-IT300	WP	Ø5.5, H=2, C=6			
	MN-IT340	NP	Ø4, H=4, C=6	Impression coping set for closed tray, internal hex.	Titanium	  
	MD-IT340	SP	Ø4.8, H=4, C=6			
	MW-IT340	WP	Ø5.5, H=4, C=6			
	MN-PF330	NP	Ø4.5, L=11	Press fit for closed tray, internal hex.	Titanium	  
	MD-PF375	SP	Ø5, L=11			
	MW-PF550	WP	Ø5.5, L=11			
	MN-RSM10	NP	Ø3.3, L=12	Implant analog, internal hex.	Stainless steel	  
	MD-RSM10	SP	Ø3.75, L=12			
	MW-RSM10	WP	Ø4.6, L=12			
	MN-S0160	NP	L=7.6	Prosthetic screw, 30 Ncm	Titanium	  
	MD-S0220	SP/WP				
	MD-G0213	SP/WP	L=23	Impression coping guide pin for open tray	Titanium	  
	MN-GP100	NP	L=16.5	Impression coping guide pin, internal hex.	Titanium	  

 Concave emergence profile H = Gingival height C = Crown height L = Length

3 Restoration Options


Single crown



Bridge



Overdenture





Temporary Restoration

Temporary Cylinders



	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MN-TC023 MD-TC023 MW-TC023	NP SP WP	Ø4, H=2, C=10 Ø4.8, H=2, C=10 Ø5.5, H=2, C=10	Temporary cylinder, concave, free rotation, internal hex.	Titanium	
	MN-TC033 MD-TC033 MW-TC033	NP SP WP	Ø4, H=3, C=10 Ø4.8, H=3, C=10 Ø5.5, H=3, C=10	Temporary cylinder, concave, free rotation, internal hex.	Titanium	
	MN-TCH23 MD-TCH23 MW-TCH23	NP SP WP	Ø4, H=2, C=10 Ø4.8, H=2, C=10 Ø5.5, H=2, C=10	Temporary cylinder, concave, anti-rotation, internal hex.	Titanium	
	MN-TCH33 MD-TCH33 MW-TCH33	NP SP WP	Ø4, H=3, C=10 Ø4.8, H=3, C=10 Ø5.5, H=3, C=10	Temporary cylinder, concave, anti-rotation, internal hex.	Titanium	
	MN-TC023 MD-TC023 MW-TC023	NP SP WP	Ø4, H=2, C=10 Ø4.8, H=2, C=10 Ø5.5, H=2, C=10	Temporary cylinder, concave, free rotation, internal hex.	Titanium	
	MN-TC033 MD-TC033 MW-TC033	NP SP WP	Ø4, H=3, C=10 Ø4.8, H=3, C=10 Ø5.5, H=3, C=10	Temporary cylinder, concave, free rotation, internal hex.	Titanium	
	MN-TCH23 MD-TCH23 MW-TCH23	NP SP WP	Ø4, H=2, C=10 Ø4.8, H=2, C=10 Ø5.5, H=2, C=10	Temporary cylinder, concave, anti-rotation, internal hex.	Titanium	
	MN-TCH33 MD-TCH33 MW-TCH33	NP SP WP	Ø4, H=3, C=10 Ø4.8, H=3, C=10 Ø5.5, H=3, C=10	Temporary cylinder, concave, anti-rotation, internal hex.	Titanium	

Concave emergence profile

H = Gingival height

C = Crown height

L = Length

	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MN-TC013	NP	Ø4, H=1, C=10	Temporary cylinder, free rotation, internal hex.	Titanium	  
	MD-TC013	SP	Ø4.8, H=1, C=10			
	MW-TC013	WP	Ø5.5, H=1, C=10			
	MN-TCH13	NP	Ø4, H=1, C=10	Temporary cylinder, anti-rotation, internal hex.	Titanium	  
	MD-TCH13	SP	Ø4.8, H=1, C=10			
	MW-TCH13	WP	Ø5.5, H=1, C=10			
	MN-TPA10	NP	Ø4, H=1, C=10	Temporary plastic abutment, anti-rotation, internal hex.	PEEK	  
	MD-TPA10	SP	Ø4.75, H=1, C=10			
	MW-TPA10	WP	Ø5.5, H=1, C=10			
	MN-TPA12	NP	Ø3.5, H=1, C=10	Narrow profile temporary plastic abutment, anti-rotation, internal hex.	PEEK	  
	MD-TPA12	SP	Ø3.75, H=1, C=10			
	MW-TPA12	WP	Ø4.5, H=1, C=10			
	MD-TPH50	SP	Ø5, H=2.5, C=9.5	Temporary plastic cylinder, anti-rotation, internal hex.	PEEK	 
	MW-TPH50	WP	Ø5.5, H=2.5, C=9.5			
	MN-S0160	NP	L=7.6	Prosthetic screw, 30 Ncm	Titanium	  
	MD-S0220	SP/WP				

 Concave emergence profile

H = Gingival height

C = Crown height

L = Length



2

Abutment Options

Ti

Titanium

PEEK

Polyether ether ketone



Cemented Restoration

Multi-Unit System



	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MA-NC248 MA-SC248	NP SP	Ø4.8, H=2	Multi-Unit concave abutment, internal hex.	Titanium	 
	MA-NC348 MA-SC348 MA-WC348	NP SP WP	Ø4.8, H=3	Multi-Unit concave abutment, internal hex.	Titanium	  
	MA-NC448 MA-SC448	NP SP	Ø4.8, H=4	Multi-Unit concave abutment, internal hex.	Titanium	 
	MA-WC548	WP	Ø4.8, H=5	Multi-Unit concave abutment, internal hex.	Titanium	
	MA-N1480 MA-S1480 MA-W1480	NP SP WP	Ø4.8, H=1	Multi-Unit abutment, internal hex.	Titanium	  
	MA-N2480 MA-S2480	NP SP	Ø4.8, H=2	Multi-Unit abutment, internal hex.	Titanium	 
	MA-N3480 MA-S3480 MA-W3480	NP SP WP	Ø4.8, H=3	Multi-Unit abutment, internal hex.	Titanium	  
	MA-N4480 MA-S4480	NP SP	Ø4.8, H=4	Multi-Unit abutment, internal hex.	Titanium	 

 Concave emergence profile

H = Gingival height

C = Crown height

L = Length

	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MA-N5480 MA-S5480 MA-W5480	NP SP WP	Ø4.8, H=5	Multi-Unit abutment, internal hex.	Titanium	  
	MA-NO171 MA-SO171 MA-WO171	NP SP WP	Ø4.8, H=1	Multi-Unit angulated abutment, 17°, internal hex.	Titanium	  
	MA-NO172 MA-SO172 MA-WO172	NP SP WP	Ø4.8, H=2	Multi-Unit angulated abutment, 17°, internal hex.	Titanium	  
	MA-SO173 MA-WO173	SP WP	Ø4.8, H=3	Multi-Unit angulated abutment, 17°, internal hex.	Titanium	 
	MA-NO301 MA-SO301 MA-WO301	NP SP WP	Ø4.8, H=1	Multi-Unit angulated abutment, 30°, internal hex.	Titanium	  
	MA-NO302 MA-SO302 MA-WO302	NP SP WP	Ø4.8, H=2	Multi-Unit angulated abutment, 30°, internal hex.	Titanium	  
	MA-SO303 MA-WO303	SP WP	Ø4.8, H=3	Multi-Unit angulated abutment, 30°, internal hex.	Titanium	 
	MU-BGC48		L=3.3	Gold bar coping abutment for Multi-Unit	Gold-e	
	MU-RSM48		L=12	Analog for Multi-Unit	Titanium	
	MU-CC480		Ø4.8 L=4.3	Cementing cap for Multi-Unit	Titanium	
	MU-CC481		Ø4.8 L=4.75	Cementing cap for Multi-Unit (New design)	Titanium	
	MU-CP048		L=10	Plastic cylinder for Multi-Unit	Burnout plastic	
	MU-H4480		Ø4.8 L=4.3	Healing cap for Multi-Unit	Titanium	

 Concave emergence profile

H = Gingival height

C = Crown height

L = Length

Multi-Unit System

	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MU-I0480		Ø4.8 L=10	Impression coping for Multi-Unit, open tray	Titanium	
	MU-IT480		Ø4.8 L=8	Impression coping for Multi-Unit, closed tray	Titanium	
	MU-TO480		Ø4.8 L=10	Titanium temporary cylinder for Multi-Unit	Titanium	
	MU-SP102		L=10	Multi-Unit scan post	Titanium	
	MU-MAN10		Ø4.8 L=10	Digital model analog for Multi-Unit	Titanium	
	MT-MURL2		L=21	Straight Multi-Unit long ratchet key	Stainless steel	
	MT-MURS2		L=16	Straight Multi-Unit short ratchet key	Stainless steel	
	MT-MUML2		L=31	Straight Multi-Unit long motor key	Stainless steel	
	MT-MUMS2		L=24	Straight Multi-Unit short motor key	Stainless steel	
	MT-RDL30		L=27	Long driver for 0.05 inch hex.	Stainless steel	
	MU-S0220		L=3.8	Prosthetic screw for Multi-Unit, 25Ncm	Titanium	
	MU-G0220		L=3.8	Prosthetic gold screw for Multi-Unit, 20Ncm	Gold-s	



Screw Retained Restoration

Plastic Cylinders



	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MN-CPH13	NP	Ø3.5, C=11	Plastic cylinder, anti-rotation, internal hex.	Burnout Plastic	  
	MD-CPH13	SP	Ø3.75, C=11			
	MW-CPH13	WP	Ø4.5, C=11			
	MN-CP013	NP	Ø3.5, C=11	Plastic cylinder, free rotation, internal hex.	Burnout Plastic	  
	MD-CP013	SP	Ø3.75, C=11			
	MW-CP013	WP	Ø4.5, C=11			
	MD-CPH50	SP	Ø5, C=11	Wide profile plastic cylinder, anti-rotation, internal hex.	Burnout Plastic	 
	MW-CPH60	WP	Ø6.5, C=11			
	MD-CP050	SP	Ø5, C=11	Wide profile plastic cylinder, free rotation, internal hex.	Burnout Plastic	 
	MW-CP060	WP	Ø6.5, C=11			

 Concave emergence profile

H = Gingival height

C = Crown height

L = Length

Plastic Cylinders

	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MN-GPC10	NP	Ø4.3, H=1, C=14	Gold plastic cylinder, anti-rotation, internal hex.	Burnout Plastic + Gold-a	  
	MD-GPC10	SP	Ø5, H=1, C=14			
	MW-GPC10	WP	Ø5.5, H=1, C=14			
	MN-GP010	NP	Ø4.3, H=1, C=14	Gold plastic cylinder, free rotation, internal hex.	Burnout plastic + Gold-a	  
	MD-GP010	SP	Ø5, H=1, C=14			
	MW-GP010	WP	Ø5.5, H=1, C=14			
	MN-CCI10	NP	Ø4.3, H=1, C=14	CoCr plastic cylinder, anti-rotation, internal hex.	Cobalt-Chrome	  
	MD-CCI10	SP	Ø5, H=1, C=14			
	MW-CCI10	WP	Ø5.5, H=1, C=14			
	MN-CC010	NP	Ø4.3, H=1, C=14	CoCr plastic cylinder, free rotation, internal hex.	Cobalt-Chrome	  
	MD-CC010	SP	Ø5, H=1, C=14			
	MW-CC010	WP	Ø5.5, H=1, C=14			
	MN-S0160	NP	L=7.6	Prosthetic screw, 30 Ncm	Titanium	  
	MD-S0220	SP/WP				
	MD-S0224	SP/WP	L=14	Long prosthetic screw, 30 Ncm	Titanium	 
	MD-S0222	SP/WP	L=12	Medium prosthetic screw, 30 Ncm	Titanium	 
	MD-G0220	SP/WP	L=7.7	Gold prosthetic screw, 20 Ncm	Gold-s	 

 Concave emergence profile H = Gingival height C = Crown height L = Length

3 Platforms

 **Narrow Platform**


 **Standard Platform**

 **Wide Platform**



Digital Restoration

EZ-Base System



Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
 MN-EZB15 MD-EZB15 MW-EZB15	NP SP WP	Ø4.25, H=1.5, C=5 Ø4.3, H=1.5, C=5 Ø5.7, H=1.5, C=5	EZ-Base abutment, internal hex.	Titanium	  
	MT-ELM10	L=30	Long motor key for angled screw channel	Stainless steel	
	MT-ESM10	L=24	Short motor key for angled screw channel	Stainless steel	
	MT-ELR10	L=22.5	Long ratchet key for angled screw channel	Stainless steel	
	MT-ESR10	L=18.7	Short ratchet key for angled screw channel	Stainless steel	
	MN-E0160 MD-E0220	L=7.1 L=7.6	EZ-Base prosthetic screw	Titanium	 

 Concave emergence profile

H = Gingival height

C = Crown height

L = Length



Digital Restoration

Ti-Bases



	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MN-TBC15 MD-TBC15 MW-TBC15	NP SP WP	Ø4.2, H=1.5, C=4 Ø4.3, H=1.5, C=4 Ø5.7, H=1.5, C=4	Ti-Base, concave, anti-rotation, internal hex.	Titanium	
	MN-FTB15 MD-FTB15 MW-FTB15	NP SP WP	Ø4.2, H=1.5, C=4 Ø4.3, H=1.5, C=4 Ø5.7, H=1.5, C=4	Ti-Base, concave, free rotation, internal hex.	Titanium	
	MN-TBC30 MD-TBC30 MW-TBC30	NP SP WP	Ø4.2, H=3, C=4 Ø4.3, H=3, C=4 Ø5.7, H=3, C=4	Ti-Base, concave, anti-rotation, internal hex.	Titanium	
	MN-FTB30 MD-FTB30 MW-FTB30	NP SP WP	Ø4.2, H=3, C=4 Ø4.3, H=3, C=4 Ø5.7, H=3, C=4	Ti-Base, concave, free rotation, internal hex.	Titanium	
	MN-TB001 MD-TB001 MW-TB001	NP SP WP	Ø4.2, H=0.5, C=4 Ø4.3, H=0.5, C=4 Ø5.7, H=0.5, C=4	Ti-Base, anti-rotation, internal hex.	Titanium	
	MN-TB002 MD-TB002 MW-TB002	NP SP WP	Ø4, H=0.5, C=4 Ø4.3, H=0.5, C=4 Ø5.7, H=0.5, C=4	Ti-Base, free rotation, internal hex.	Titanium	
	MN-IN060 MD-IN060	NP SP	Ø4.2, H=0.5, C=6 Ø4.3, H=0.5, C=6	Incisor esthetic Ti abutment, internal hex.	Titanium	
	MN-S0160 MD-S0220	NP SP/WP	L=7.6	Prosthetic screw, 30 Ncm	Titanium	

Concave emergence profile

H = Gingival height

C = Crown height

L = Length



Digital Restoration

Scan Posts



MK-0056
Scan post kit, internal hex.

	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MN-SP102 MD-SP102 MW-SP102	NP SP WP	L=10	Scan post, internal hex.	Titanium	  
	MU-SP102		L=10	Multi-Unit scan post	Titanium	
	MN-S0160 MD-S0220	NP SP/WP	L=7.6	Prosthetic screw, 20 Ncm	Titanium	  

 Concave emergence profile

H = Gingival height


C = Crown height

L = Length



Digital Restoration

Digital Model Analogs



	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MN-MAN10	NP	Ø3.30, L=10	Digital model analog, internal hex.	Titanium	  
	MD-MAN10	SP	Ø3.80, L=10			
	MW-MAN10	WP	Ø4.60, L=10			
	MU-MAN10		Ø4.80 L=10	Digital model analog for Multi-Unit	Titanium	

 Concave emergence profile

H = Gingival height

C = Crown height

L = Length



Digital Restoration

Titanium Blanks



	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MN-IB091 MD-IB091 MW-IB091	NP SP WP	Ø9	Titanium Blank, anti-rotation, internal hex.	Titanium	  
	MN-IB121 MD-IB121 MW-IB121	NP SP WP	Ø12	Titanium Blank, anti-rotation, internal hex.	Titanium	  
	MN-AGT12 MD-AGT12 MW-AGT12	NP SP WP	Ø12	Titanium Blank for Amann Girrbach, anti-rotation, internal hex.	Titanium	  
	MN-MAB10 MD-MAB10 MW-MAB10	NP SP WP	Ø5	Milling analog for Titanium Blank, internal hex.	Stainless steel	  

 Concave emergence profile

H = Gingival height


C = Crown height

L = Length



Cemented Restoration

Cementable Abutments



Cat. No. Platform Dimensions (mm) Description Material Platform



MN-AT102

NP

Ø4, H=2, C=7

Cementable angulated abutment, 10°, concave, internal hex.

Titanium

NP



MD-AT152
MW-AT152

SP
WP

Ø4.8, H=2, C=7.5
Ø5.5, H=2, C=7.5

Cementable angulated abutment, 15°, concave, internal hex.

Titanium

SP WP



MD-AT153

SP

Ø4.8, H=3, C=7.5

Cementable angulated abutment, 15°, concave, internal hex.

Titanium

SP



MN-AT202

NP

Ø4, H=2, C=7

Cementable angulated abutment, 20°, concave, internal hex.

Titanium

NP

Concave emergence profile

H = Gingival height

C = Crown height

L = Length

	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MD-AT252	SP	Ø4.8, H=2, C=7.5	Cementable angulated abutment, 25°, concave, internal hex.	Titanium	
	MD-AT253	SP	Ø4.8, H=3, C=7.5	Cementable angulated abutment, 25°, concave, internal hex.	Titanium	
	MD-A0010 MW-P0010	SP WP	Ø5.5, H=1, C=8 Ø6.3, H=1, C=8	Cementable abutment, internal hex.	Titanium	 
	MD-P0030	SP	Ø5.6, H=3, C=8	Cementable abutment, internal hex.	Titanium	
	MN-AT101	NP	Ø4, H=1, C=7	Cementable angulated abutment, 10°, internal hex.	Titanium	
	MD-AN150	SP	Ø4.8, H=1, C=6.4	Cementable angulated abutment, 15°, internal hex.	Titanium	
	MD-A1510 MW-P1510	SP WP	Ø5.2, H=1, C=8 Ø5.8, H=1, C=8	Cementable angulated abutment, 15°, internal hex.	Titanium	 
	MN-AN020	NP	Ø4, H=1, C=7.3	Cementable angulated abutment, 20°, internal hex.	Titanium	
	MD-AN151 MW-AN151	SP WP	Ø4.9, H=1, C=7.5 Ø6, H=1, C=7.3	Cementable angulated abutment, 15°, internal hex.	Titanium	 
	MD-P1530	SP	Ø5.2, H=3, C=8	Cementable angulated abutment, 15°, internal hex.	Titanium	
	MD-A2510	SP	Ø5.6, H=1, C=8	Cementable angulated abutment, 25°, internal hex.	Titanium	
	MD-AN251	SP	Ø5.9, H=1, C=7	Cementable angulated abutment, 25°, internal hex.	Titanium	

 Concave emergence profile

H = Gingival height

C = Crown height


L = Length

Cementable
Abutments

	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MD-P2530	SP	Ø5.7, H=3, C=8	Cementable angulated abutment, 25°, internal hex.	Titanium	
	MN-S0160 MD-S0220	NP SP/WP	L=7.6	Prosthetic screw, 30 Ncm	Titanium	  
	MD-S0224	SP/WP	L=14	Long prosthetic screw, 30 Ncm	Titanium	 
	MD-G0220	SP/WP	L=7.7	Gold prosthetic screw, 20 Ncm	Gold-s	 

 Concave emergence profile H = Gingival height C = Crown height L = Length

 5 Angle Options





Cemented Restoration

Cementing Posts



	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MN-MC102	NP	Ø4, H=2, C=10	Cementing post, concave, internal hex.	Titanium	  
	MD-MC102	SP	Ø4.8, H=2, C=10			
	MW-MC102	WP	Ø5.5, H=2, C=10			
	MN-MAC10	NP	Ø4, H=1, C=10	Cementing post, internal hex.	Titanium	  
	MD-MAC10	SP	Ø4.75, H=1, C=10			
	MW-MAC10	WP	Ø5.5, H=1, C=10			
	MD-WMAC1	SP	Ø5.5, H=1, C=10	Wide profile cementing post, internal hex.	Titanium	 
	MW-WMAC1	WP	Ø6.5, H=1, C=10			
	MD-MACF1	SP	Ø4.8, H=2.5, C=8.5	Friction fit cementing post with platform switching, internal hex.	Titanium	
	MD-CTP10	SP	Ø3.75, H=0, C=8.5	Narrow profile cementing post, internal hex.	Titanium	 
	MW-CTP10	WP	Ø5.5, H=0, C=8.5			
	MN-S0160	NP	L=7.6	Prosthetic screw, 30 Ncm	Titanium	  
	MD-S0220	SP/WP				
	MD-S0224	SP/WP	L=14	Long prosthetic screw, 30 Ncm	Titanium	 
	MD-S0222	SP/WP	L=12	Medium prosthetic screw, 30 Ncm	Titanium	 
	MD-G0220	SP/WP	L=7.7	Gold prosthetic screw, 20 Ncm	Gold-s	 

 Concave emergence profile

H = Gingival height

C = Crown height

L = Length



Cemented Restoration

CPK - Complete
Prosthetic Kit

	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MK-NPC62	NP	Ø4, H=2, C=6	Complete Prosthetic Kit, concave, internal hex.	Titanium	
	MK-CPC62	SP	Ø4.8, H=2, C=6			
	MK-WPC62	WP	Ø5.5, H=2, C=6			
	MK-NPC63	NP	Ø4, H=3, C=6	Complete Prosthetic Kit, concave, internal hex.	Titanium	
	MK-CPC63	SP	Ø4.8, H=3, C=6			
	MK-WPC63	WP	Ø5.5, H=3, C=6			
	MK-NPC64	NP	Ø4, H=4, C=6	Complete Prosthetic Kit, concave, internal hex.	Titanium	
	MK-CPC64	SP	Ø4.8, H=4, C=6			
	MK-WPC64	WP	Ø5.5, H=4, C=6			
	MK-NPK41	NP	Ø4, H=1, C=4	Complete Prosthetic Kit, internal hex.	Titanium	
	MK-CPK41	SP	Ø4.8, H=1, C=4			
	MK-WPK41	WP	Ø5.5, H=1, C=4			
	MK-NPK42	NP	Ø4, H=2, C=4	Complete Prosthetic Kit, internal hex.	Titanium	
	MK-CPK42	SP	Ø4.8, H=2, C=4			
	MK-WPK42	WP	Ø5.5, H=2, C=4			
	MK-NPK43	NP	Ø4, H=3, C=4	Complete Prosthetic Kit, internal hex.	Titanium	
	MK-CPK43	SP	Ø4.8, H=3, C=4			
	MK-WPK43	WP	Ø5.5, H=3, C=4			
	MK-NPK44	NP	Ø4, H=4, C=4	Complete Prosthetic Kit, internal hex.	Titanium	
	MK-CPK44	SP	Ø4.8, H=4, C=4			
	MK-WPK44	WP	Ø5.5, H=4, C=4			
	MK-NPK61	NP	Ø4, H=1, C=6	Complete Prosthetic Kit, internal hex.	Titanium	
	MK-CPK61	SP	Ø4.8, H=1, C=6			
	MK-WPK61	WP	Ø5.5, H=1, C=6			
	MK-NPK62	NP	Ø4, H=2, C=6	Complete Prosthetic Kit, internal hex.	Titanium	
	MK-CPK62	SP	Ø4.8, H=2, C=6			
	MK-WPK62	WP	Ø5.5, H=2, C=6			

Concave emergence profile

H = Gingival height

C = Crown height

L = Length

	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MK-NPK63	NP	Ø4, H=3, C=6	Complete Prosthetic Kit, internal hex.	Titanium	  
	MK-CPK63	SP	Ø4.8, H=3, C=6			
	MK-WPK63	WP	Ø5.5, H=3, C=6			
	MK-NPK64	NP	Ø4, H=4, C=6	Complete Prosthetic Kit, internal hex.	Titanium	  
	MK-CPK64	SP	Ø4.8, H=4, C=6			
	MK-WPK64	WP	Ø5.5, H=4, C=6			
	MK-CPK81	SP	Ø4.8, H=1, C=8	Complete Prosthetic Kit, internal hex.	Titanium	 
	MK-WPK81	WP	Ø5.5, H=1, C=8			
	MK-CPK82	SP	Ø4.8, H=2, C=8	Complete Prosthetic Kit, internal hex.	Titanium	 
	MK-WPK82	WP	Ø5.5, H=2, C=8			
	MK-CPK83	SP	Ø4.8, H=3, C=8	Complete Prosthetic Kit, internal hex.	Titanium	 
	MK-WPK83	WP	Ø5.5, H=3, C=8			
	MK-CPK84	SP	Ø4.8, H=4, C=8	Complete Prosthetic Kit, internal hex.	Titanium	 
	MK-WPK84	WP	Ø5.5, H=4, C=8			
	MN-CPC62	NP	Ø4, H=2, C=6	Transgingival abutment, concave, internal hex.	Titanium	  
	MD-CPC62	SP	Ø4.8, H=2, C=6			
	MW-CPC62	WP	Ø5.5, H=2, C=6			
	MN-CPC63	NP	Ø4, H=3, C=6	Transgingival abutment, concave, internal hex.	Titanium	  
	MD-CPC63	SP	Ø4.8, H=3, C=6			
	MW-CPC63	WP	Ø5.5, H=3, C=6			
	MN-CPC64	NP	Ø4, H=4, C=6	Transgingival abutment, concave, internal hex.	Titanium	  
	MD-CPC64	SP	Ø4.8, H=4, C=6			
	MW-CPC64	WP	Ø5.5, H=4, C=6			
	MN-CPK41	NP	Ø4, H=1, C=4	Transgingival abutment, internal hex.	Titanium	  
	MD-CPK41	SP	Ø4.8, H=1, C=4			
	MW-CPK41	WP	Ø5.5, H=1, C=4			
	MN-CPK42	NP	Ø4, H=2, C=4	Transgingival abutment, internal hex.	Titanium	  
	MD-CPK42	SP	Ø4.8, H=2, C=4			
	MW-CPK42	WP	Ø5.5, H=2, C=4			
	MN-CPK43	NP	Ø4, H=3, C=4	Transgingival abutment, internal hex.	Titanium	  
	MD-CPK43	SP	Ø4.8, H=3, C=4			
	MW-CPK43	WP	Ø5.5, H=3, C=4			

 Concave emergence profile

H = Gingival height

C = Crown height

L = Length

CPK - Complete Prosthetic Kit

	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MN-CPK44	NP	Ø4, H=4, C=4	Transgingival abutment, internal hex.	Titanium	  
	MD-CPK44	SP	Ø4.8, H=4, C=4			
	MW-CPK44	WP	Ø5.5, H=4, C=4			
	MN-CPK61	NP	Ø4, H=1, C=6	Transgingival abutment, internal hex.	Titanium	  
	MD-CPK61	SP	Ø4.8, H=1, C=6			
	MW-CPK61	WP	Ø5.5, H=1, C=6			
	MN-CPK62	NP	Ø4, H=2, C=6	Transgingival abutment, internal hex.	Titanium	  
	MD-CPK62	SP	Ø4.8, H=2, C=6			
	MW-CPK62	WP	Ø5.5, H=2, C=6			
	MN-CPK63	NP	Ø4, H=3, C=6	Transgingival abutment, internal hex.	Titanium	  
	MD-CPK63	SP	Ø4.8, H=3, C=6			
	MW-CPK63	WP	Ø5.5, H=3, C=6			
	MN-CPK64	NP	Ø4, H=4, C=6	Transgingival abutment, internal hex.	Titanium	  
	MD-CPK64	SP	Ø4.8, H=4, C=6			
	MW-CPK64	WP	Ø5.5, H=4, C=6			
	MD-CPK81	SP	Ø4.8, H=1, C=8	Transgingival abutment, internal hex.	Titanium	 
	MW-CPK81	WP	Ø5.5, H=1, C=8			
	MD-CPK82	SP	Ø4.8, H=2, C=8	Transgingival abutment, internal hex.	Titanium	 
	MW-CPK82	WP	Ø5.5, H=2, C=8			
	MD-CPK83	SP	Ø4.8, H=3, C=8	Transgingival abutment, internal hex.	Titanium	 
	MW-CPK83	WP	Ø5.5, H=3, C=8			
	MD-CPK84	SP	Ø4.8, H=4, C=8	Transgingival abutment, internal hex.	Titanium	 
	MW-CPK84	WP	Ø5.5, H=4, C=8			
	MN-IC600	NP	Ø4.2, L=8	Transfer plastic cap for CPK	Plastic	  
	MD-IC800	SP	Ø5, L=10			
	MW-IC800	WP	Ø5.7, L=10			
	MM-CHC40	NP	C=4	Plastic healing cap	PEEK	  
	MM-CHC60	SP	C=6			
	MM-CHC80	WP	C=8			
	MN-RSM60	NP	Ø4, C=6, L=18	Abutment analog for CPK	Stainless steel	  
	MD-RSM60	SP	Ø4.8, C=6, L=18			
	MW-RSM60	WP	Ø5.5, C=6, L=18			



Overdentures

LOCKiT Attachment System



	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MN-LOC01	NP	Ø3.75, H=1	LOCKiT abutment, internal hex.	Titanium alloy + Titanium nitride	  
	MD-LOC01	SP	Ø3.75, H=1			
	MW-LOC01	WP	Ø4.5, H=1			
	MN-LOC02	NP	Ø3.75, H=2	LOCKiT abutment, internal hex.	Titanium alloy + Titanium nitride	 
	MD-LOC02	SP	Ø3.75, H=2			
	MN-LOC03	NP	Ø3.75, H=3	LOCKiT abutment, internal hex.	Titanium alloy + Titanium nitride	  
	MD-LOC03	SP	Ø3.75, H=3			
	MW-LOC03	WP	Ø4.5, H=3			
	MN-LOC04	NP	Ø3.75, H=4	LOCKiT abutment, internal hex.	Titanium alloy + Titanium nitride	 
	MD-LOC04	SP	Ø3.75, H=4			
	MN-LOC05	NP	Ø3.75, H=5	LOCKiT abutment, internal hex.	Titanium alloy + Titanium nitride	  
	MD-LOC05	SP	Ø3.75, H=5			
	MW-LOC05	WP	Ø4.5, H=5			

 Concave emergence profile

H = Gingival height

C = Crown height

L = Length

LOCKiT Attachment System

	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MK-NLC1	NP	Ø3.75, H=1	LOCKiT abutment kit, internal hex.	Titanium alloy + Titanium nitride	  
	MK-SLC1	SP	Ø3.75, H=1			
	MK-WLC1	WP	Ø4.5, H=1			
	MK-NLC2	NP	Ø3.75, H=2	LOCKiT abutment kit, internal hex.	Titanium alloy + Titanium nitride	 
	MK-SLC2	SP	Ø3.75, H=2			
	MK-NLC3	NP	Ø3.75, H=3	LOCKiT abutment kit, internal hex.	Titanium alloy + Titanium nitride	  
	MK-SLC3	SP	Ø3.75, H=3			
	MK-WLC3	WP	Ø4.5, H=3			
	MK-NLC4	NP	Ø3.75, H=4	LOCKiT abutment kit, internal hex.	Titanium alloy + Titanium nitride	 
	MK-SLC4	SP	Ø3.75, H=4			
	MK-NLC5	NP	Ø3.75, H=5	LOCKiT abutment kit, internal hex.	Titanium alloy + Titanium nitride	  
	MK-SLC5	SP	Ø3.75, H=5			
	MK-WLC5	WP	Ø4.5, H=5			
	MT-LOC10		Ø7.8 L=86	LOCKiT core tool	Stainless steel	
	MM-LOC1		Ø5.4 L=4.4	LOCKiT impression coping	Titanium	
	MM-LOA40	NP/SP	Ø4, L=10.2	LOCKiT analog	Titanium	  
	MM-LOA50	WP	Ø5, L=10.2			
	MM-LOC02			LOCKiT male processing package		
	MM-CLE10		Ø4.75 H=1.78	LOCKiT male retention cap 10° standard (clear)	Nylon	
	MM-PIN10		Ø4.72 H=1.78	LOCKiT male retention cap 10° light (pink)	Nylon	
	MM-BLU10		Ø4.72 H=1.78	LOCKiT male retention cap 10° extra light (blue)	Nylon	



Concave emergence profile

H = Gingival height

C = Crown height

L = Length

	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MM-GRE10		Ø4.75 H=1.78	LOCKiT male retention cap 20° standard (green)	Nylon	
	MM-RED10		Ø4.72 H=1.78	LOCKiT male retention cap 20° extra light (red)	Nylon	
	MM-BLA10		Ø4.75 H=1.9	LOCKiT male processing cap for impression coping (black)	Polyethylene	
	MM-TIT10		Ø5.5 H=2.25	LOCKiT metal housing titanium cap	Titanium	
	MT-RSD19		L=19	LOCKiT short torque ratchet driver	Stainless steel	
	MT-RSD24		L=24	LOCKiT long torque ratchet driver	Stainless steel	
	MT-LOC11		Ø7.35 H=4.1	LOCKiT abutment holder sleeve		



Concave emergence profile

H = Gingival height


C = Crown height

L = Length




5

Retention Options




Extra Light Retention




Light Retention



Standard Retention



Extra Light Retention



Standard Retention



Overdentures

OT-Equator Kits



	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MK-NOE1 MK-SOE1 MK-WOE1	NP SP WP	Ø4, H=1 Ø4, H=1 Ø5, H=1	OT-Equator kit, internal hex.	Titanium	NP SP WP
	MK-NOE2 MK-SOE2	NP SP	Ø4, H=2 Ø4, H=2	OT-Equator kit, internal hex.	Titanium	NP SP
	MK-NOE3 MK-SOE3 MK-WOE3	NP SP WP	Ø4, H=3 Ø4, H=3 Ø5, H=3	OT-Equator kit, internal hex.	Titanium	NP SP WP
	MK-NOE4 MK-SOE4	NP SP	Ø4, H=4 Ø4, H=4	OT-Equator kit, internal hex.	Titanium	NP SP

Concave emergence profile

H = Gingival height

C = Crown height

L = Length

	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MK-NOE5	NP	Ø4, H=5	OT-Equator kit, internal hex.	Titanium	  
	MK-SOE5	SP	Ø4, H=5			
	ET-IT002			OT-Equator handling tool		
	OE-IC002			OT-Equator impression coping cap		
	OE-RSM60			OT-Equator analog	Stainless steel	
	MT-RDL30		L=27	Long driver for 0.05 inch hex.	Stainless steel	
	MT-RDS30		L=21.3	Short driver for 0.05 inch hex.	Stainless steel	
	OE-RCP01			OT-Equator soft retentive cap (pink), 4-pack		
	OE-RCW01			OT-Equator standard retentive cap (white), 4-pack		
	OE-RCV01			OT-Equator rigid retentive cap (violet), 4-pack		
	OE-RCY01			OT-Equator extra soft retentive cap (yellow), 4-pack		
	OE-RCB01			OT-Equator laboratory retentive cap (black), 4-pack		
	OE-MH001			OT-Equator metal housing		

 Concave emergence profile

H = Gingival height


C = Crown height

L = Length



Overdentures

Ball Attachment System



	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MB-N1375 MW-B1500	SP WP	Ø4, H=1 Ø5, H=1	Ball attachment anchor, internal hex.	Titanium alloy + Titanium nitride	 
	MB-N2330 MB-N2375	NP SP	Ø4, H=2	Ball attachment anchor, internal hex.	Titanium alloy + Titanium nitride	 
	MB-N3330 MB-N3375 MW-B3500	NP SP WP	Ø4, H=3 Ø4, H=3 Ø5, H=3	Ball attachment anchor, internal hex.	Titanium alloy + Titanium nitride	  
	MB-N4330 MB-N4375	NP SP	Ø4, H=4	Ball attachment anchor, internal hex.	Titanium alloy + Titanium nitride	 
	MB-N5375 MW-B5500	SP WP	Ø4, H=5 Ø5, H=5	Ball attachment anchor, internal hex.	Titanium alloy + Titanium nitride	 
	MB-AN152	SP	H=2	Ball anchor angulated abutment, 15°, internal hex.	Titanium	
	MB-AN252	SP	H=2	Ball anchor angulated abutment, 25°, internal hex.	Titanium	

 Concave emergence profile

H = Gingival height

C = Crown height

L = Length

	Cat. No.	Platform	Dimensions (mm)	Description	Material	Platform
	MB-RS375		L=15.1	Ball attachment analog	Stainless steel	
	MB-SF200		Ø3.7 L=3.5	Plastic ball cap for ball attachment	Plastic	
	MB-GPS10		L=2.4	Gold ball cap for ball attachment	Gold-e	
	MB-TTP10		Ø3.5 L=3.1	Titanium ball cap for ball attachment	Titanium	
	MB-PR010		Ø3.8 L=2	Plastic mounting sleeve for ball attachment	Plastic	
	MB-TBS10		Ø2.5	Titanium ball cap spring for ball attachment, standard	Steel	
	MB-TBS20		Ø2.6	Titanium ball cap spring for ball attachment, soft	Steel	
	MB-DMH10			Double metal housing for ball attachment		
	MB-PPC10		Ø4.2 L=2.6	Standard plastic cap for ball attachment metal housing	Plastic	
	MB-PPR10		Ø4.2 L=2.6	Soft plastic cap for ball attachment metal housing, 800-950g retention	Plastic	
	MB-DB235		Ø15 L=2.6	Plastic disc for ball attachment	Plastic	
	MT-RB225		L=22	Ball attachment anchor key	Stainless steel	

 Concave emergence profile

H = Gingival height

C = Crown height

L = Length



MGUIDE - DIGITAL DENTISTRY

MGUIDE is an advanced yet simple solution for guided surgery. The 3D-printed template is designed with an open frame that allows visibility, irrigation and accessibility from all angles without the need for removal.

Designed as a keyless system, MGUIDE allows for single-handed procedures, eliminating the need for unnecessary tools. With MGUIDE, clinicians can choose between flapless and raised-flap procedures.





MGUIDE

MGUIDE Templates



MGUIDE prosthetic template

MGUIDE surgical template

	Cat. No	Description
	MG-TMP10	MGUIDE surgical template
	MG-TMP20	MGUIDE surgical template for edentulous cases
	MG-TMP50	MGUIDE surgical and prosthetic templates

Please note: MGUIDE prosthetic template might not be registered or available in every country/region.



MGUIDE

MGUIDE SEVEN Surgical Set



MG-K001
MGUIDE SEVEN Surgical Set

	Cat. No	Dimensions (mm)	Description	Material	Sold
	MG-TP047	Ø4.7 L=26	MGUIDE tissue punch, standard platform	Stainless steel	 
	MG-TP055	Ø5.5 L=26	MGUIDE tissue punch, wide platform	Stainless steel	 
	MG-BM047	Ø4.7 L=24.5	MGUIDE bone mill, standard platform	Stainless steel	 
	MG-BM055	Ø5.5 L=24.5	MGUIDE bone mill, wide platform	Stainless steel	 
	MG-D0624	Ø2.4 L=30.8	MGUIDE 6mm pilot drill	Stainless steel	 
	MG-D0824	Ø2.4 L=32.8	MGUIDE pilot drill for 8mm length implant	Stainless steel	 

L = Length

MGUIDE SEVEN Surgical Set

	Cat. No	Dimensions (mm)	Description	Material	Sold
	MG-D1024	Ø2.4 L=34.8	MGUIDE pilot drill for 10mm length implant	Stainless steel	 
	MG-D1124	Ø2.4 L=36.3	MGUIDE pilot drill for 11.5mm length implant	Stainless steel	 
	MG-D1324	Ø2.4 L=37.8	MGUIDE pilot drill for 13mm length implant	Stainless steel	 
	MG-D0633	Ø2.2-3.2 L=30.8	MGUIDE drill for Ø3.30 x 6mm length implant	Stainless steel	 
	MG-D0833	Ø2.2-3.2 L=32.8	MGUIDE drill for Ø3.30 x 8mm length implant	Stainless steel	 
	MG-D1033	Ø2.2-3.2 L=34.8	MGUIDE drill for Ø3.30 x 10mm length implant	Stainless steel	 
	MG-D1133	Ø2.2-3.2 L=36.3	MGUIDE drill for Ø3.30 x 11.5mm length implant	Stainless steel	 
	MG-D1333	Ø2.2-3.2 L=37.8	MGUIDE drill for Ø3.30 x 13mm length implant	Stainless steel	 
	MG-D0637	Ø2.8-3.6 L=30.8	MGUIDE drill for Ø3.75 x 6mm length implant	Stainless steel	 
	MG-D0837	Ø2.8-3.6 L=32.8	MGUIDE drill for Ø3.75 x 8mm length implant	Stainless steel	 
	MG-D1037	Ø2.8-3.6 L=34.8	MGUIDE drill for Ø3.75 x 10mm length implant	Stainless steel	 
	MG-D1137	Ø2.8-3.6 L=36.3	MGUIDE drill for Ø3.75 x 11.5mm length implant	Stainless steel	 

L = Length



Cat. No	Dimensions (mm)	Description	Material	Sold
MG-D1337	Ø2.8-3.6 L=37.8	MGUIDE drill for Ø3.75 x 13mm length implant	Stainless steel	 
MG-D0642	Ø3.3-4.1 L=30.8	MGUIDE drill for Ø4.20 x 6mm length implant	Stainless steel	 
MG-D0842	Ø3.3-4.1 L=32.8	MGUIDE drill for Ø4.20 x 8mm length implant	Stainless steel	 
MG-D1042	Ø3.3-4.1 L=34.8	MGUIDE drill for Ø4.20 x 10mm length implant	Stainless steel	 
MG-D1142	Ø3.3-4.1 L=36.3	MGUIDE drill for Ø4.20 x 11.5mm length implant	Stainless steel	 
MG-D1342	Ø3.3-4.1 L=37.8	MGUIDE drill for Ø4.20 x 13mm length implant	Stainless steel	 
MG-D0650	Ø4.1-4.9 L=30.8	MGUIDE drill for Ø5 x 6mm length implant	Stainless steel	 
MG-D0850	Ø4.1-4.9 L=32.8	MGUIDE drill for Ø5 x 8mm length implant	Stainless steel	 
MG-D1050	Ø4.1-4.9 L=34.8	MGUIDE drill for Ø5 x 10mm length implant	Stainless steel	 
MG-D1150	Ø4.1-4.9 L=36.3	MGUIDE drill for Ø5 x 11.5mm length implant	Stainless steel	 
MG-D1350	Ø4.1-4.9 L=37.8	MGUIDE drill for Ø5 x 13mm length implant	Stainless steel	 
MG-DLG55	Ø5.5 L=55	MGUIDE drill length gauge	Titanium	 

L = Length

MGUIDE SEVEN
Surgical Set

	Cat. No	Dimensions (mm)	Description	Material	Sold
	MG-GRI01	Ø7 L=29.7	MGUIDE direct ratchet insertion tool, internal hex.	Stainless steel	 
	MG-GRN01	Ø7 L=29.7	MGUIDE direct ratchet insertion tool, internal hex., narrow platform	Stainless steel	 
	MG-GMI10	Ø5.4 L=28.9	MGUIDE motor insertion tool, internal hex.	Stainless steel	 
	MG-GRI10	Ø7 L=24.5	MGUIDE ratchet insertion tool, internal hex.	Stainless steel	 
	MT-RDL30	L=27	Long driver for 0.05 inch hex.	Stainless steel	 
	MG-IE172	L=32.7	MGUIDE insertion tool extractor, internal hex., standard/wide platform	Titanium	 
	MG-IE160	L=32.7	MGUIDE insertion tool extractor, internal hex., narrow platform	Titanium	 
	MG-DFP20	Ø2 L=35.3	MGUIDE drill for fixation pin	Stainless steel	 
	MG-FP020	Ø2 L=27	MGUIDE fixation pin	Titanium	 
	MG-TAS55	Ø5.5 L=24	MGUIDE template anchoring screw	Titanium	 
	MT-RA480		Ratchet to hex. adapter	Stainless steel	 
	MT-RI030	L=84	Monoblock ratchet		 

L = Length



MGUIDE

MGUIDE Drill Replacement Kit

MG-SD01 
MGUIDE Drill Replacement Kit



Cat. No	Dimensions (mm)	Description	Material	Sold
 MG-TP047	Ø4.7 L=26	MGUIDE tissue punch, standard platform	Stainless steel	 
 MG-TP055	Ø5.5 L=26	MGUIDE tissue punch	Stainless steel	 
 MG-BM047	Ø4.7 L=24.5	MGUIDE bone mill, standard platform	Stainless steel	 
 MG-BM055	Ø5.5 L=24.5	MGUIDE bone mill, wide platform	Stainless steel	 
 MG-D0624	Ø2.4 L=30.8	MGUIDE pilot drill for 6mm length implant	Stainless steel	 
 MG-D0824	Ø2.4 L=32.8	MGUIDE pilot drill for 8mm length implant	Stainless steel	 
 MG-D1024	Ø2.4 L=34.8	MGUIDE pilot drill for 10mm length implant	Stainless steel	 
 MG-D1124	Ø2.4 L=36.3	MGUIDE pilot drill for 11.5mm length implant	Stainless steel	 
 MG-D1324	Ø2.4 L=37.8	MGUIDE pilot drill for 13mm length implant	Stainless steel	 

L = Length



MGUIDE

MGUIDE SEVEN
Surgical Set, Narrow
Sleeve

MG-K002
MGUIDE SEVEN Surgical Set,
Narrow Sleeve



	Cat. No	Dimensions (mm)	Description	Material	Sold
  	MG-NBM40	Ø4 L=24.5	MGUIDE Narrow Sleeve bone mill	Stainless steel	 
	MG-NTP40	Ø4 L=26	MGUIDE Narrow Sleeve tissue punch	Stainless steel	 
	MG-N0624	Ø2.4 L=30.8	MGUIDE Narrow Sleeve pilot drill for 6mm length implant	Stainless steel	 

L = Length

Cat. No	Dimensions (mm)	Description	Material	Sold
 MG-N0824	Ø2.4 L=32.8	MGUIDE Narrow Sleeve pilot drill for 8mm length implant	Stainless steel	 
 MG-N1024	Ø2.4 L=34.8	MGUIDE Narrow Sleeve pilot drill for 10mm length implant	Stainless steel	 
 MG-N1124	Ø2.4 L=36.3	MGUIDE Narrow Sleeve pilot drill for 11.5mm length implant	Stainless steel	 
 MG-N1324	Ø2.4 L=37.8	MGUIDE Narrow Sleeve pilot drill for 13mm length implant	Stainless steel	 
 MG-N0633	Ø2.2-3.2 L=30.8	MGUIDE Narrow Sleeve drill for Ø3.30 x 6mm length implant	Stainless steel	 
 MG-N0833	Ø2.2-3.2 L=32.8	MGUIDE Narrow Sleeve drill for Ø3.30 x 8mm length implant	Stainless steel	 
 MG-N1033	Ø2.2-3.2 L=34.8	MGUIDE Narrow Sleeve drill for Ø3.30 x 10mm length implant	Stainless steel	 
 MG-N1133	Ø2.2-3.2 L=36.3	MGUIDE Narrow Sleeve drill for Ø3.30 x 11.5mm length implant	Stainless steel	 
 MG-N1333	Ø2.2-3.2 L=37.8	MGUIDE Narrow Sleeve drill for Ø3.30 x 13mm length implant	Stainless steel	 
 MG-N0637	Ø2.8-3.6 L=30.8	MGUIDE Narrow Sleeve drill for Ø3.75 x 6mm length implant	Stainless steel	 
 MG-N0837	Ø2.8-3.6 L=32.8	MGUIDE Narrow Sleeve drill for Ø3.75 x 8mm length implant	Stainless steel	 
 MG-N1037	Ø2.8-3.6 L=34.8	MGUIDE Narrow Sleeve drill for Ø3.75 x 10mm length implant	Stainless steel	 

L = Length

MGUIDE SEVEN
Surgical Set, Narrow
Sleeve

	Cat. No	Dimensions (mm)	Description	Material	Sold
	MG-N1137	Ø2.8-3.6 L=36.3	MGUIDE Narrow Sleeve drill for Ø3.75 x 11.5mm length implant	Stainless steel	 
	MG-N1337	Ø2.8-3.6 L=37.8	MGUIDE Narrow Sleeve drill for Ø3.75 x 13mm length implant	Stainless steel	 
	MG-NLG40	Ø4 L=55	MGUIDE Narrow Sleeve drill length gauge	Titanium	 
	MG-NMN10	Ø5 L=30.4	MGUIDE Narrow Sleeve motor insertion tool, internal hex., narrow platform	Stainless steel	 
	MG-NMS10	Ø5 L=28.9	MGUIDE Narrow Sleeve motor insertion tool, internal hex., standard platform	Stainless steel	 
	MG-NRN01	Ø5.3 L=29.5	MGUIDE Narrow Sleeve direct ratchet insertion tool, internal hex., narrow platform	Stainless steel	 
	MG-NRS01	Ø5.3 L=24.7	MGUIDE Narrow Sleeve direct ratchet insertion tool, internal hex., standard platform	Stainless steel	 
	MG-NRN10	Ø5.3 L=25	MGUIDE Narrow Sleeve ratchet insertion tool, internal hex., narrow platform	Stainless steel	 
	MG-NRS10	Ø5.3 L=24.5	MGUIDE Narrow Sleeve ratchet insertion tool, internal hex., standard platform	Stainless steel	 
	MG-NTAS0	Ø5.3 L=24	MGUIDE Narrow Sleeve template anchoring screw	Titanium	 

L = Length



MGUIDE

MGUIDE Drill Replacement Kit, Narrow Sleeve

MG-ND01 
MGUIDE Drill Replacement Kit,
Narrow Sleeve



Cat. No	Dimensions (mm)	Description	Material	Sold
 MG-NBM40	Ø4 L=24.5	MGUIDE Narrow Sleeve bone mill	Stainless steel	 
 MG-NTP40	Ø4 L=26	MGUIDE Narrow Sleeve tissue punch	Stainless steel	 
 MG-N0624	Ø2.4 L=30.8	MGUIDE Narrow Sleeve pilot drill for 6mm length implant	Stainless steel	 
 MG-N0824	Ø2.4 L=32.8	MGUIDE Narrow Sleeve pilot drill for 8mm length implant	Stainless steel	 
 MG-N1024	Ø2.4 L=34.8	MGUIDE Narrow Sleeve pilot drill for 10mm length implant	Stainless steel	 
 MG-N1124	Ø2.4 L=36.3	MGUIDE Narrow Sleeve pilot drill for 11.5mm length implant	Stainless steel	 
 MG-N1324	Ø2.4 L=37.8	MGUIDE Narrow Sleeve pilot drill for 13mm length implant	Stainless steel	 

L = Length




ADDITIONAL INFORMATION



Additional information


Pilot Drills





Additional information


SEVEN Twist Drills





Additional information


M4 Twist Drills





Additional information

SEVEN and M4 Short Twist Drills





Additional information

Alloy Composition

	Material	Titanium Alloy Ti 6Al 4V ELI	Gold-a	Gold-s	Gold-e
	Application	C, K	L, G, K	L, G, K	L, K
Composition	Gold %	-	60	61	70
	Platinum %	-	19	13.5	8.50
	Silver %	-	-	16.5	13.40
	Palladium %	-	20	-	-
	Copper %	-	-	9	7.50
	Zinc %		-	Balance	0.05
	Iridium %		Balance	-	0.10
	Titanium (max)	Balance	-	-	-
	Carbon (max)	0.08	-	-	-
	Iron (max)	0.25	-	-	-
	Oxygen (max)	0.13	-	-	-
	Nitrogen (max)	0.05	-	-	-
	Hydrogen (max)	0.013	-	-	-
	Aluminum	5.5-6.5	-	-	-
	Vanadium	3.5-4.5	-	-	-
	Melting range °C	1604-1660	1400-1490	950-1050	895-1010
Hardness	Soft annealing HV5	-	145	160	170
	Hardened HV5	350-385	240	250	295
	Soft °C / min	-	100/30	750/5	750/5
	Hardened °C / min	-	700/30	350/15	450/15

Legend: C–bonding K–resin L–soldering G–casting onto





Additional information

INDEX

Cat. No. Page

C

CT-BTC24	45
CT-BTC24	48
CT-P2406	30
CT-P2406	47
CT-P2408	30
CT-P2408	47
CT-P2408	51
CT-P2410	30
CT-P2410	47
CT-P2410	51
CT-P2411	31
CT-P2411	47
CT-P2411	51
CT-P2413	31
CT-P2413	48
CT-P2413	51
CT-P2416	31
CT-P2416	45
CT-P2416	48
CT-P2416	51
CT-P2416	52
CT-UN001	61
CT-US001	61
CT-UW001	62

Cat. No. Page

E

ET-IT002	44
ET-IT002	62
ET-IT002	91

M

MA-N1480	70
MA-N2480	70
MA-N3480	70
MA-N4480	70
MA-N5480	71
MA-NC248	70
MA-NC348	70
MA-NC448	70
MA-NO171	71
MA-NO172	71
MA-NO301	71
MA-NO302	71
MA-S1480	70
MA-S2480	70
MA-S3480	70
MA-S4480	70
MA-S5480	71
MA-SC248	70
MA-SC348	70

Cat. No. Page

MA-SC448	70
MA-SO171	71
MA-SO172	71
MA-SO173	71
MA-SO301	71
MA-SO302	71
MA-SO303	71
MA-W1480	70
MA-W3480	70
MA-W5480	71
MA-WC348	70
MA-WC548	70
MA-WO171	71
MA-WO172	71
MA-WO173	71
MA-WO301	71
MA-WO302	71
MA-WO303	71
MB-AN152	92
MB-AN252	92
MB-DB235	93
MB-DMH10	93
MB-GPS10	93
MB-N1375	92
MB-N2330	92
MB-N2375	92
MB-N3330	92
MB-N3375	92
MB-N4330	92
MB-N4375	92
MB-N5375	92
MB-PPC10	93
MB-PPR10	93
MB-PR010	93
MB-RS375	93
MB-SF200	93
MB-TBS10	93
MB-TBS20	93

Cat. No.	Page
MB-TTP10	93
MD-A0010	81
MD-A1510	81
MD-A2510	81
MD-AGT12	79
MD-AN150	81
MD-AN151	81
MD-AN251	81
MD-AT152	80
MD-AT153	80
MD-AT252	81
MD-AT253	81
MD-CC010	74
MD-CCI10	74
MD-CP013	73
MD-CP050	73
MD-CPC62	85
MD-CPC63	85
MD-CPC64	85
MD-CPH13	73
MD-CPH50	73
MD-CPK41	85
MD-CPK42	85
MD-CPK43	85
MD-CPK44	86
MD-CPK61	86
MD-CPK62	86
MD-CPK63	86
MD-CPK64	86
MD-CPK81	86
MD-CPK82	86
MD-CPK83	86
MD-CPK84	86
MD-CTP10	83
MD-E0220	75
MD-EZB15	75
MD-FTB15	76
MD-FTB30	76
MD-G0213	67
MD-G0220	74
MD-G0220	82
MD-G0220	83
MD-GP010	74
MD-GPC10	74
MD-I0375	66
MD-IB091	79
MD-IB121	79
MD-IC800	86
MD-IN060	76
MD-IO481	66
MD-IT100	66

Cat. No.	Page
MD-IT300	67
MD-IT340	67
MD-ITC30	66
MD-LOC01	87
MD-LOC02	87
MD-LOC03	87
MD-LOC04	87
MD-LOC05	87
MD-MAB10	79
MD-MAC10	83
MD-MACF1	83
MD-MAN10	78
MD-MC102	83
MD-P0030	81
MD-P1530	81
MD-P2530	82
MD-PF375	50
MD-PF375	67
MD-RSM10	67
MD-RSM60	86
MD-S0220	67
MD-S0220	69
MD-S0220	74
MD-S0220	76
MD-S0220	77
MD-S0220	82
MD-S0220	83
MD-S0222	74
MD-S0222	82
MD-S0222	83
MD-S0224	74
MD-S0224	82
MD-S0224	83
MD-SP102	77
MD-TB001	76
MD-TB002	76
MD-TBC15	76
MD-TBC30	76
MD-TC013	69
MD-TC023	68
MD-TC023	68
MD-TC033	68
MD-TC033	68
MD-TCH13	69
MD-TCH23	68
MD-TCH23	68
MD-TCH33	68
MD-TCH33	68
MD-TPA10	69
MD-TPA12	69
MD-TPH50	69

Cat. No.	Page
MD-WMAC1	83
MF4-06420	23
MF4-06500	25
MF4-06600	27
MF4-08375	21
MF4-08420	23
MF4-08500	25
MF4-08600	27
MF4-10330	19
MF4-10375	21
MF4-10420	23
MF4-10500	25
MF4-10600	27
MF4-11330	19
MF4-11375	21
MF4-11420	23
MF4-11500	25
MF4-11600	27
MF4-13330	19
MF4-13375	21
MF4-13420	23
MF4-13500	25
MF4-13600	27
MF4-16330	19
MF4-16375	21
MF4-16420	23
MF4-16500	25
MF7-06600	15
MF7-08600	15
MF7-10600	15
MF7-11600	15
MF7-13600	15
MF7-D06420	11
MF7-D06500	13
MF7-D08375	9
MF7-D08420	11
MF7-D08500	13
MF7-D10330	7
MF7-D10375	9
MF7-D10420	11
MF7-D10500	13
MF7-D11330	7
MF7-D11375	9
MF7-D11420	11
MF7-D11500	13
MF7-D13330	7
MF7-D13375	9
MF7-D13420	11
MF7-D13500	13
MF7-D16330	7
MF7-D16375	9

Cat. No.	Page
MF7-D16420	11
MF7-D16500	13
MG1-TE480	44
MG-BM047	97
MG-BM047	101
MG-BM055	97
MG-BM055	101
MG-D0624	97
MG-D0624	101
MG-D0633	98
MG-D0637	98
MG-D0642	99
MG-D0650	99
MG-D0824	97
MG-D0824	101
MG-D0833	98
MG-D0837	98
MG-D0842	99
MG-D0850	99
MG-D1024	98
MG-D1024	101
MG-D1033	98
MG-D1037	98
MG-D1042	99
MG-D1050	99
MG-D1124	98
MG-D1124	101
MG-D1133	98
MG-D1137	98
MG-D1142	99
MG-D1150	99
MG-D1324	98
MG-D1324	101
MG-D1333	98
MG-D1337	99
MG-D1342	99
MG-D1350	99
MG-DFP20	100
MG-DLG55	99
MG-FP020	100
MG-GMI10	100
MG-GRI01	100
MG-GRI10	100
MG-GRN01	100
MG-IE160	100
MG-IE172	100
MG-K001	97
MG-K002	102
MG-N0624	102
MG-N0624	105
MG-N0633	103

Cat. No.	Page
MG-N0637	103
MG-N0824	103
MG-N0824	105
MG-N0833	103
MG-N0837	103
MG-N1024	103
MG-N1024	105
MG-N1033	103
MG-N1037	103
MG-N1124	103
MG-N1124	105
MG-N1133	103
MG-N1137	104
MG-N1324	103
MG-N1324	105
MG-N1333	103
MG-N1337	104
MG-NBM40	102
MG-NBM40	105
MG-ND01	105
MG-NLG40	104
MG-NMN10	104
MG-NMS10	104
MG-NRN01	104
MG-NRN10	104
MG-NRS01	104
MG-NRS10	104
MG-NTAS0	104
MG-NTP40	102
MG-NTP40	105
MG-SD01	101
MG-TAS55	100
MG-TMP10	96
MG-TMP20	96
MG-TMP50	96
MG-TP047	97
MG-TP047	101
MG-TP055	97
MG-TP055	101
MK-0001	63
MK-0016	52
MK-0020	58
MK-0039	60
MK-0056	77
MK-0059	55
MK-BS001	56
MK-CPC62	84
MK-CPC63	84
MK-CPC64	84
MK-CPK41	84
MK-CPK42	84

Cat. No.	Page
MK-CPK43	84
MK-CPK44	84
MK-CPK61	84
MK-CPK62	84
MK-CPK63	85
MK-CPK64	85
MK-CPK81	85
MK-CPK82	85
MK-CPK83	85
MK-CPK84	85
MK-NLC1	88
MK-NLC2	88
MK-NLC3	88
MK-NLC4	88
MK-NLC5	88
MK-NOE1	90
MK-NOE2	90
MK-NOE3	90
MK-NOE4	90
MK-NOE5	91
MK-NPC62	84
MK-NPC63	84
MK-NPC64	84
MK-NPK41	84
MK-NPK42	84
MK-NPK43	84
MK-NPK44	84
MK-NPK61	84
MK-NPK62	84
MK-NPK63	85
MK-NPK64	85
MK-SD01	51
MK-SDS06	57
MK-SDS(06-13)	57
MK-SDS08	57
MK-SDS10	57
MK-SDS11	57
MK-SDS13	57
MK-SLC1	88
MK-SLC2	88
MK-SLC3	88
MK-SLC4	88
MK-SLC5	88
MK-SOE1	90
MK-SOE2	90
MK-SOE3	90
MK-SOE4	90
MK-SOE5	91
MK-SVXD	45
MK-T048	47
MK-WLC1	88

Cat. No. Page

MK-WLC3	88
MK-WLC5	88
MK-WOE1	90
MK-WOE3	90
MK-WPC62	84
MK-WPC63	84
MK-WPC64	84
MK-WPK41	84
MK-WPK42	84
MK-WPK43	84
MK-WPK44	84
MK-WPK61	84
MK-WPK62	84
MK-WPK63	85
MK-WPK64	85
MK-WPK81	85
MK-WPK82	85
MK-WPK83	85
MK-WPK84	85
MK-XD01	45
MM-BLA10	89
MM-BLU10	88
MM-CHC40	86
MM-CHC60	86
MM-CHC80	86
MM-CLE10	88
MM-GRE10	89
MM-LOA40	88
MM-LOA50	88
MM-LOC02	88
MM-LOC11	88
MM-PIN10	88
MM-RED10	89
MM-TIT10	89
MN-AGT12	79
MN-AN020	81
MN-AT101	81
MN-AT102	80
MN-AT202	80
MN-CC010	74
MN-CCI10	74
MN-CP013	73
MN-CPC62	85
MN-CPC63	85
MN-CPC64	85
MN-CPH13	73
MN-CPK41	85
MN-CPK42	85
MN-CPK43	85
MN-CPK44	86
MN-CPK61	86

Cat. No. Page

MN-CPK62	86
MN-CPK63	86
MN-CPK64	86
MN-E0160	75
MN-EZB15	75
MN-FTB15	76
MN-FTB30	76
MN-GP010	74
MN-GP100	67
MN-GPC10	74
MN-I0330	66
MN-IB091	79
MN-IB121	79
MN-IC600	86
MN-IN060	76
MN-IO401	66
MN-IT100	66
MN-IT300	67
MN-IT340	67
MN-ITC30	66
MN-LOC01	87
MN-LOC02	87
MN-LOC03	87
MN-LOC04	87
MN-LOC05	87
MN-MAB10	79
MN-MAC10	83
MN-MAN10	78
MN-MC102	83
MN-PF330	50
MN-PF330	67
MN-RSM10	67
MN-RSM60	86
MN-S0160	67
MN-S0160	69
MN-S0160	74
MN-S0160	76
MN-S0160	77
MN-S0160	82
MN-S0160	83
MN-SP102	77
MN-TB001	76
MN-TB002	76
MN-TBC15	76
MN-TBC30	76
MN-TC013	69
MN-TC023	68
MN-TC023	68
MN-TC033	68
MN-TC033	68
MN-TCH13	69

Cat. No. Page

MN-TCH23	68
MN-TCH23	68
MN-TCH33	68
MN-TCH33	68
MN-TPA10	69
MN-TPA12	69
MO-CO160	58
MO-CO200	58
MO-CO240	58
MO-CO280	58
MO-SO300	58
MO-SO350	58
MO-SO400	58
MT-00400	38
MT-00500	38
MT-00600	38
MT-BGP32	55
MT-BGP37	55
MT-BGP46	55
MT-BP320	55
MT-BP370	55
MT-BP460	55
MT-BTI20	40
MT-BTT28	48
MT-BTT32	48
MT-BTT40	48
MT-BTT45	48
MT-BTT50	49
MT-CLR21	60
MT-CS330	37
MT-CS375	37
MT-CS420	37
MT-CS500	37
MT-CSN33	37
MT-CSN33	46
MT-CSN33	53
MT-DE001	41
MT-DE001	46
MT-DE001	49
MT-DE001	53
MT-ELM10	44
MT-ELM10	75
MT-ELR10	44
MT-ELR10	62
MT-ELR10	75
MT-ESM10	44
MT-ESM10	75
MT-ESR10	44
MT-ESR10	75
MT-GDN33	37
MT-GDN33	46

Cat. No. Page

MT-GDN33	49
MT-GDN33	53
MT-GDN50	37
MT-GDN50	46
MT-GDN50	49
MT-GDN50	53
MT-HLI10	39
MT-HLI10	46
MT-HLI10	49
MT-HLI10	53
MT-HLI10	59
MT-HLI21	39
MT-HLI21	46
MT-HLI21	49
MT-HSI10	39
MT-HSI10	49
MT-HSI10	53
MT-HSI10	59
MT-HSI21	39
MT-HTS10	43
MT-IT100	42
MT-LM005	41
MT-LM005	42
MT-LM005	53
MT-LM005	60
MT-LOC10	44
MT-LOC10	62
MT-LOC10	88
MT-LOC11	89
MT-LRH20	39
MT-LRH20	46
MT-LRH20	49
MT-LRH20	54
MT-LRH20	59
MT-LRH21	39
MT-LRH21	46
MT-LRH21	49
MT-LRH21	53
MT-MFXGS	41
MT-MSD20	42
MT-MSD20	62
MT-MUML2	43
MT-MUML2	61
MT-MUML2	72
MT-MUMS2	44
MT-MUMS2	72
MT-MURL2	43
MT-MURL2	61
MT-MURL2	72
MT-MURS2	43
MT-MURS2	72

Cat. No. Page

MT-PD440	36
MT-PDM24	45
MT-PP240	41
MT-PP240	49
MT-PP240	53
MT-RA480	41
MT-RA480	44
MT-RA480	100
MT-RB225	43
MT-RB225	61
MT-RB225	93
MT-RDE30	40
MT-RDE30	43
MT-RDL30	40
MT-RDL30	43
MT-RDL30	46
MT-RDL30	50
MT-RDL30	54
MT-RDL30	60
MT-RDL30	72
MT-RDL30	91
MT-RDL30	100
MT-RDM30	41
MT-RDM30	43
MT-RDS30	40
MT-RDS30	43
MT-RDS30	46
MT-RDS30	50
MT-RDS30	60
MT-RDS30	91
MT-RE160	61
MT-RE172	61
MT-RI030	41
MT-RI030	54
MT-RI030	59
MT-RI030	100
MT-RI040	44
MT-RI040	62
MT-RMR10	46
MT-RMR10	49
MT-RMR10	54
MT-RSD19	44
MT-RSD19	61
MT-RSD19	89
MT-RSD24	44
MT-RSD24	61
MT-RSD24	89
MT-RT070	41
MT-RT070	46
MT-RT070	50
MT-SM005	41

Cat. No. Page

MT-SM005	42
MT-SM005	60
MT-SMD10	36
MT-SMD10	47
MT-SMD10	52
MT-SMD10	59
MT-SRA10	41
MT-SRH20	39
MT-SRH20	49
MT-SRH20	54
MT-SRH20	59
MT-SRH21	39
MT-TDI28	34
MT-TDI32	34
MT-TDI38	34
MT-TDI40	34
MT-TDI45	35
MT-TDI50	35
MT-TDI55	35
MT-TDN19	36
MT-TDN19	45
MT-TDN19	47
MT-TDN19	52
MT-TDN28	33
MT-TDN28	52
MT-TDN32	33
MT-TDN32	52
MT-TDN38	33
MT-TDN38	53
MT-TDN45	33
MT-TDN45	53
MT-TDN50	33
MT-TDN50	53
MT-TDN55	33
MT-TDN55	53
MT-TDT28	32
MT-TDT28	48
MT-TDT28	51
MT-TDT32	32
MT-TDT32	48
MT-TDT32	51
MT-TDT40	32
MT-TDT40	48
MT-TDT40	51
MT-TDT45	32
MT-TDT45	48
MT-TDT45	51
MT-TDT50	32
MT-TDT50	48
MT-TP420	41
MT-TP500	41

Cat. No. Page

MT-UN002	43
MT-UN002	61
MT-US002	43
MT-US002	61
MT-UW002	43
MT-UW002	61
MU-BGC48	71
MU-CC480	71
MU-CC481	71
MU-CP048	71
MU-G0220	72
MU-H4480	71
MU-I0480	72
MU-IT480	72
MU-MAN10	72
MU-MAN10	78
MU-RSM48	71
MU-S0220	72
MU-SP102	72
MU-SP102	77
MU-TO480	72
MW-AGT12	79
MW-AN151	81
MW-AT152	80
MW-B1500	92
MW-B3500	92
MW-B5500	92
MW-CC010	74
MW-CCI10	74
MW-CP013	73
MW-CP060	73
MW-CPC62	85
MW-CPC63	85
MW-CPC64	85
MW-CPH13	73
MW-CPH60	73
MW-CPK41	85
MW-CPK42	85
MW-CPK43	85
MW-CPK44	86
MW-CPK61	86
MW-CPK62	86
MW-CPK63	86
MW-CPK64	86
MW-CPK81	86
MW-CPK82	86
MW-CPK83	86
MW-CPK84	86
MW-CTP10	83
MW-EZB15	75
MW-FTB15	76
MW-FTB30	76

Cat. No. Page

MW-GP010	74
MW-GPC10	74
MW-I0470	66
MW-IB091	79
MW-IB121	79
MW-IC800	86
MW-IO551	66
MW-IT100	66
MW-IT300	67
MW-IT340	67
MW-ITC30	66
MW-LOC01	87
MW-LOC03	87
MW-LOC05	87
MW-MAB10	79
MW-MAC10	83
MW-MAN10	78
MW-MC102	83
MW-P0010	81
MW-P1510	81
MW-PF550	67
MW-RSM10	67
MW-RSM60	86
MW-SP102	77
MW-TB001	76
MW-TB002	76
MW-TBC15	76
MW-TBC30	76
MW-TC013	69
MW-TC023	68
MW-TC023	68
MW-TC033	68
MW-TC033	68
MW-TCH13	69
MW-TCH23	68
MW-TCH23	68
MW-TCH33	68
MW-TCH33	68
MW-TPA10	69
MW-TPA12	69
MW-TPH50	69
MW-WMAC1	83

Cat. No. Page

OE-IC002	91
OE-MH001	91
OE-RCB01	91
OE-RCP01	91
OE-RCV01	91
OE-RCW01	91
OE-RCY01	91
OE-RSM60	91



