





Conical Connection Prosthetic Options

MIS prosthetic components are suitable for a wide range of clinical applications, including implant level or abutment level, one-stage or two-stage treatments, and single or multiple-unit solutions. They come in various gingival heights, diameters, and angles.



Additional information may be added to the product descriptions without prior notice.
© MIS Implants Technologies Ltd. All rights reserved.



NOT FOR USE IN THE USA AND CANADA

Products may not be available in some regulated markets. Please contact your local Dentsply Sirona sales office for current product assortment and availability.

The MIS Quality System complies with international quality standard ISO 13485: 2016 - Quality Management System for Medical Devices, with Medical Device Directive 93/42/EEC, and with EU Medical Device Regulation MDR 2017/745. Please note, that not all products are registered or available in every country/region.

IFUs for MIS products may be found at: <https://ifu.mis-implants.com>. Adobe Acrobat is required to view the IFU file on the website. This software may be freely downloaded from the Adobe website.

Table of Contents

Implants
Dimensions

P. 6 - 9

Prosthetic Options
Listed by Healing Caps

P. 10 - 15

Healing Caps

P. 16 - 18

Impression
Coping

P. 19 - 20

Temporary
Restoration

P. 21 - 23

○ Temporary Cylinders

Screw-retained Restoration

P. 24 - 34

- CONNECT System
- Multi-Unit System
- Gold Plastic Cylinders
- EZ-Base System

Digital Restoration

P. 35 - 40

- Ti-Base
- Digital Model Analogs
- Scan Posts
- Titanium Blanks

Cemented Restoration

P. 41 - 49

- Cementable Abutments
- CPK

Overdenture

P. 50 - 55

- LOCKIT System
- OT-Equator Kits
- Ball Attachment System



Prosthetic Screws and Drivers

P. 56 - 57



Implant Dimensions: MIS LYNX

Please note, MIS LYNX might
not be registered or available
in every country/region.





Narrow
Platform



Standard
Platform




Wide
Platform





Implant Dimensions: C1



Narrow
Platform



Standard
Platform






Wide
Platform







Implant Dimensions: V3

Narrow
Platform



Standard
Platform





NARROW PLATFORM

Prosthetic Options Listed by Healing Caps


MIS LYNX / C1/ V3 Implants

Each column lists the different heights of healing caps and corresponding prosthetic options. **Bold** item numbers indicate the products shown in the images above.


Item numbers beginning with **VN** are for **V3** implants only.

Item numbers beginning with **CN** are for **MIS LYNX** and **C1** implants.


1. Thin healing cap



2. Standard healing cap




3. Anatomic healing cap



1.

Healing
caps

Impression
coping,
open tray



Gingival heights 3, 4, 5, 6, 8 mm

CN-H0333

CN-H0433

CN-H0533

CN-H0633

CN-H0833

CN-IO335
CN-IO405
CN-IO331
CN-IO330
VN-IO401
VN-IO331
VN-IO335

CN-IO335
CN-IO405


2.

Healing caps	Impression coping, open tray	Impression coping, closed tray	Cemented options	Screw retained options
--------------	------------------------------	--------------------------------	------------------	------------------------

or

or

or



Gingival heights 3, 4, 5, 6, 8 mm	CN-HS340	CN-IO335 VN-IO405 CN-IO331 CN-IO330 VN-IO401 VN-IO405 VN-IO331 VN-IO335	CN-IC401 VN-IC401 CN-IC405 VN-IC405 CN-PF330 VN-PF340 CN-IT300 VN-IC400	CN-MAC10 VN-MC101 CK-NPK41 VK-NPK41 CK-NPK61 VK-NPK61	CK-NOE1 CN-B1330
	CN-HS440			CK-NPK42 VK-NPK42 CK-NPK62 VK-NPK62	CK-NOE2 CN-B2330
	CN-HS540		CN-IC401 VN-IC401 CN-IC405 VN-IC405 CN-PF330 VN-PF340	CK-NPK43 VK-NPK43 CK-NPK63 VK-NPK63 CN-MAC30 VN-MC103	CK-NOE3 CN-B3330
	CN-HS640	CN-IO335 CN-IO405 VN-IO335 VN-IO405		CK-NPK44 VK-NPK44 CK-NPK64 VK-NPK64	CK-NOE4 CN-B4330
	CN-HS840				CK-NOE5


3.

Healing caps	Impression coping, open tray	Impression coping, closed tray	Cemented options	Screw retained options	Digital options
--------------	------------------------------	--------------------------------	------------------	------------------------	-----------------

or

or

or



Gingival heights 2, 3, 4, 5, 6, 8 mm	CN-HA248	CN-IO335 VN-IO405 CN-IO331 CN-IO330 VN-IO401 VN-IO405 VN-IO331 VN-IO335	CN-IC401 VN-IC401 CN-IC405 VN-IC405 CN-PF330 VN-PF340 CN-IT300 VN-IC400		CN-TB001 VN-TB001 CN-TB002 CN-IN060 VN-IN060
	CN-HA348			CN-AN001 VN-AT001 CN-AN010 VN-AT101 CN-AN020 VN-AT201	CN-GP010 CN-GPC10 VN-GPC10 CM-N1480
	CN-HA448			CN-AT002 VN-AT002 CN-AT102 VN-AT102 CN-AT202 VN-AT202	CM-N2480
	CN-HA548	CN-IO335 CN-IO405 CN-IO331 CN-IO330 VN-IO335 VN-IO405 VN-IO331 VN-IO401	CN-IC401 VN-IC401 CN-IC405 VN-IC405 CN-PF330 VN-PF340	CN-AN003 VN-AT003 CN-AN310 VN-AT103 CN-AN320 VN-AT203	CN-TB030 VN-TB030 CN-TBF30 CN-IN360 VN-IN360
	CN-HA648	CN-IO335 VN-IO335 CN-IO405 VN-IO405			CM-N4480
	CN-HA848				CM-N5480


STANDARD PLATFORM

Prosthetic Options Listed by Healing Caps


MIS LYNX / C1/ V3 Implants

Each column lists the different heights of healing caps and corresponding prosthetic options. **Bold** item numbers indicate the products shown in the images above.


1. Thin healing cap



2. Standard healing cap




3. Anatomic healing cap



1.

Healing
caps

Impression
coping,
open tray



Gingival heights 3, 4, 5, 6, 8 mm

CS-H0339

CS-H0439

CS-H0539

CS-H0639

CS-H0839

CS-IO395
CS-IO391

CS-IO395

2.

Healing caps

Impression coping, open tray


Impression coping, closed tray

Cemented options

Screw retained options

Digital options

Screw retained options



Gingival heights 3, 4, 5, 6, 8 mm

CS-HS248

CS-IO375
CS-IO485
CS-IO391
CS-IO395CS-IC481
CS-IC485
CS-PF375
CS-IT300CS-TB001
CS-IN060
CS-TB002

CS-HS348

CS-A0010
CS-A1510
CS-A2510
CS-MAC10
CK-CPK41
CK-CPK61
CK-CPK81CM-S1480
CS-GPC10
CS-GP010CS-TB015
CS-IN160
CS-TBF15
CS-EZB15
CS-LEZ15CK-SOE1
CS-B1375

CS-HS448

CS-AT002
CS-AT152
CS-AT252
CK-CPK42
CK-CPK62
CK-CPK82

CM-S2480

CK-SOE2
CS-B2375**CS-HS548**CS-IC481
CS-IC485
CS-PF375CS-MC103
CS-A0030
CS-A1530
CS-A2530
CK-CPK43
CK-CPK63
CK-CPK83**CM-S3480****CS-TB030**
CS-TBF30
CS-IN360**CK-SOE3**
CS-B3375

CS-HS648

CS-IO485
CS-IO395CS-IC485
CS-PF375CK-CPK44
CK-CPK64
CK-CPK84

CM-S4480

CK-SOE4
CS-B4375

CS-HS848

CM-S5480

CK-SOE5
CS-B5375


3.

Healing caps

Impression coping, open tray

Impression coping, closed tray

Cemented options



Gingival heights 2, 3, 4, 5, 6, 8 mm

VS-HS358

CS-IO375
CS-IO485
CS-IO391
CS-IO395CS-IC481
CS-IC485
CS-PF375
CS-IT300VS-AT001
CS-AN151
CS-AN251
VS-MC101
VK-CPK61
VK-CPK81**VS-HS458**CS-IC481
CS-IC485
CS-IT300
CS-PF375**VS-AT002**
VS-AT152
VS-AT252
VK-CPK62
VK-CPK82

VS-HS558

VS-AT003
VS-AT153
VS-AT253
VS-MC103
VK-CPK63
VK-CPK83

VS-HS658

CS-MC103
CS-A0030
CS-A1530
CS-A2530
CK-CPK43
CK-CPK63
CK-CPK83

VS-HS858


CS-IO485
CS-IO395VK-CPK64
VK-CPK84

CS-HS848

2.

Healing caps	Impression coping, open tray	Impression coping, closed tray	Cemented options
--------------	------------------------------	--------------------------------	------------------

or



Gingival heights 2, 3, 5 mm	CW-HA263	CW-IO500 CW-IO555 CW-IO501 CW-IO505	CW-IC551 CW-IC555 CW-PF500 CW-IT300	CW-TB001 CW-TB002
	CW-HA363			CW-TB015 CW-TBF15
	CW-HA563		CW-IC551 CW-IC555 CW-PF500	CW-TB030 CW-TBF30

Healing Caps Narrow Platform

Tightening torque: 20-25 Ncm

	HEIGHT						
		MIS LYNX/ C1/V3 2mm	MIS LYNX/ C1/V3 3mm	MIS LYNX/ C1/V3 4mm	MIS LYNX/ C1/V3 5mm	MIS LYNX/ C1/V3 6mm	MIS LYNX/ C1/V3 8mm
Thin, Ø3.3							
Standard, Ø4							
Anatomic, Ø4.8							

Healing Caps Standard Platform

Tightening torque: 20-25 Ncm

HEIGHT

Thin, Ø3.9

Standard, Ø4.8

Anatomic, Ø5.8

MIS LYNX/
C1/V3

2mm

MIS LYNX/
C1/V3

3mm

MIS LYNX/
C1/V3

4mm

MIS LYNX/
C1/V3

5mm

MIS LYNX/
C1/V3

6mm

MIS LYNX/
C1/V3

7/8mm



CS-H0339



CS-H0439



CS-H0539



CS-H0639



CS-H0839



CS-HS248




CS-HS348



CS-HS448



CS-HS548




CS-HS648



CS-HS848




VS-HS358



VS-HS458



VS-HS558



VS-HS658



VS-HS858

Healing Caps Wide Platform

Tightening torque: 20-25 Ncm

		MIS LYNX/ C1 2mm	MIS LYNX/ C1 3mm	MIS LYNX/ C1 4mm	MIS LYNX/ C1 5mm	MIS LYNX/ C1 6mm
HEIGHT						
Standard, Ø5.5						
			CW-HS355	CW-HS455	CW-HS555	CW-HS655
	Anatomic, Ø6.3					
		CW-HA263	CW-HA363	CW-HA563		

Open tray

PLATFORM:

Impression coping, 11 mm

Impression coping, 15mm

Thin impression coping, 11 mm

Thin impression coping, 15mm

Implant analog, open/closed tray



















MIS LYNX/
C1
NP

V3

NP

MIS LYNX/
C1/V3
SP

MIS LYNX/
C1
WP






Impression Coping

Closed tray



Tightening torque: 20-25 Ncm

PLATFORM:



Impression coping, 11mm




Impression coping, 15mm






Impression coping set, 2mm



Impression coping set, 4mm



Press fit



Temporary Restoration





Temporary Cylinders

Free rotation





Tightening torque: 20-25 Ncm

PLATFORM:





Temporary cylinder, 1 mm



Temporary cylinder, 2mm



Temporary cylinder, 3mm



MIS LYNX/
C1/V3
NP

MIS LYNX/
C1/V3
SP

MIS LYNX/
C1/V3
SP

MIS LYNX/
C1
WP

Temporary Restoration






Temporary Cylinders

Anti rotation






Tightening torque: 20-25 Ncm

PLATFORM






Temporary cylinder, 1mm



Temporary cylinder, 2mm



Temporary cylinder, 3mm



MIS LYNX/
C1
NP

V3
NP

MIS LYNX/
C1/V3
SP

MIS LYNX/
C1/V3
SP

MIS LYNX/
C1
WP

Temporary Restoration

Temporary Plastic Cylinders

Anti rotation

Tightening torque: 20-25 Ncm

PLATFORM

Temporary plastic abutment, 1.5mm





Temporary plastic cylinder

MIS LYNX/
C1
NP

V3
NP

MIS LYNX/
C1/V3
SP

MIS LYNX/
C1
WP



Screw Retained Restoration

CONNECT System

Ø4mm

Tightening torque: 30 Ncm

PLATFORM.

i CONNECT products shown in this catalog are sold subject to individual regulatory certification only.

CONNECT abutment, 1.5mm

CONNECT abutment, 2mm

CONNECT abutment, 3mm

CONNECT abutment, 4mm

MIS LYNX/
C1/V3

NP



CN-C2040



CN-C3040

MIS LYNX/
C1/V3


SP



CS-C1540



CS-C2040




CS-C3040




CS-C4040

MIS LYNX/
C1


WP




CW-C1540



CW-C2040



CW-C3040



CW-C4040

CONNECT System Components


Ø4mm

Healing caps



MM-H0540
MM-H1540
MM-H3040


Abutment analog



MM-RSM41

Impression coping for
closed tray,
free/anti-rotation

Requires MT-IT100 key




○ MM-IC040
○ MM-ICI40

Long impression coping
for open tray,
free/anti-rotation,



○ MM-IO011
○ MM-IOI11

Short impression coping
for open tray,
free/anti-rotation



○ MM-IO040
○ MM-IOI40

Temporary cylinders,
free/anti-rotation



○ MM-TC041
○ MM-TCI41

Final esthetic abutments,
free/anti-rotation,
height 6mm




○ MM-CE046
○ MM-CEI46

Final esthetic abutments,
free/anti-rotation,
height 8mm



○ MM-CF048
○ MM-CFI48

CoCr plastic cylinders,
free/anti-rotation



○ MM-CCR04
○ MM-CCRI4

Angulated
cementable abutment,
20°,anti-rotation



MM-AN204

Scan post,
anti-rotation



MM-SP104

Digital model abalog



MM-MAN40

Motor insertion tool



MT-CLM21 - Long
MT-CSM21 - Short

Ratchet insertion tool



MT-CLR21 - Long
MT-CSR21 - Short

Screw Retained Restoration

CONNECT System

Ø5.7mm

Tightening torque: 30 Ncm

i CONNECT products shown in this catalog are sold subject to individual regulatory certification only.

PLATFORM:

CONNECT abutment, 1.5mm

CONNECT abutment, 2mm


CONNECT abutment, 3mm

CONNECT abutment, 4mm

MIS LYNX/
C1/V3
SP



CS-C1557



CS-C2057



CS-C3057



CS-C4057

MIS LYNX/
C1
WP



CW-C1557



CW-C2057



CW-C3057



CW-C4057

CONNECT System Components Ø5.7mm

Healing caps



MM-H0557
MM-H1557
MM-H3057

Abutment analog



MM-RSM57

Impression coping
for closed tray, free/
anti-rotation

Requires MT-IT100 key



○ MM-IC057
○ MM-ICI57

Impression coping for
open tray,
free/anti-rotation



○ MM-IO057
○ MM-IOI57

Temporary cylinders,
free/anti-rotation




○ MM-TC057
○ MM-TCI57

Final esthetic abutments,
free/anti-rotation,
height 6mm



○ MM-CE056
○ MM-CEI56

Final esthetic abutments,
free/anti-rotation,
height 4mm



○ MM-CF054
○ MM-CFI54

CoCr plastic cylinders,
free/anti-rotation



○ MM-CCR05
○ MM-CCRi5

Scan post,
anti-rotation



MM-SP105

Digital model abalog



MM-MAN57

Motor insertion tool



MT-CLM21 - Long
MT-CSM21 - Short

Ratchet insertion tool



MT-CLR21 - Long
MT-CSR21 - Short

Screw Retained Restoration

Multi-Unit System

Tightening torque: 30 Ncm

PLATFORM.

Multi-Unit, 1mm

MIS LYNX/
C1/V3

NP



CM-N1480

MIS LYNX/
C1/V3

SP



CM-S1480


MIS LYNX/
C1

WP




CM-W1480

Multi-Unit, 2mm




CM-N2480



CM-S2480


Multi-Unit, 3mm



CM-N3480



CM-S3480



CM-W3480

Multi-Unit, 4mm



CM-S4480

Multi-Unit, 5mm



CM-N5480



CM-S5480

Screw Retained Restoration

Multi-Unit System

Tightening torque: 20 Ncm

PLATFORM


17° Angulated Multi-Unit, 1mm

30° Angulated Multi-Unit, 1mm

17° Angulated Multi-Unit, 2mm

30° Angulated Multi-Unit, 2mm

MIS LYNX/
C1
NP




CM-NO171



CM-NO301




CM-NO172




CM-NO302


V3
NP



VM-NO171



VM-NO301



VM-NO172



VM-NO302

Multi-Unit System


Tightening torque: 20 Ncm

PLATFORM


MIS LYNX/
C1/V3

SP

17° / 30° Angulated Multi-Unit, 1mm




CM-SO171




CM-SO301

17° / 30° Angulated Multi-Unit, 2mm




CM-SO172




CM-SO302

17° / 30° Angulated Multi-Unit, 3mm



CM-SO173



CM-SO303

Screw Retained Restoration



Multi-Unit System

Tightening torque: 20 Ncm



PLATFORM

MIS LYNX/
C1
WP



17° / 30° Angulated Multi-Unit, 1mm



17° / 30° Angulated Multi-Unit, 2mm



17° / 30° Angulated Multi-Unit, 3mm



Multi-Unit System Components

Healing cap




MU-H4480

Analog



MU-RSM48

Impression coping,
open tray



MU-I0480

Impression coping,
closed tray



MU-IT480

Temporary cylinder



MU-TO480

Plastic cylinder



MU-CP048

Gold bar coping abutment



MU-BGC48

Cementing cap
(New design)



MU-CC481

Cementing cap



MU-CC480

Scan Post



MU-SP102

Digital model analog



MU-MAN10

Screw Retained Restoration

Gold Plastic Cylinders

Free / anti rotation

Tightening torque: 30 Ncm

PLATFORM

MIS LYNX/
C1

NP

V3

NP





MIS LYNX/
C1/V3

SP



MIS LYNX/
C1

WP

Gold plastic cylinder, 1 mm,
anti-rotation



Gold plastic cylinder, 1 mm,
free rotation



Screw Retained Restoration

EZ-Base System

Tightening torque: 20-25 Ncm

PLATFORM

EZ-Base abutment, 1.5mm









EZ-Base adjustable abutment, 1.5mm

MIS LYNX/
C1
NP

V3
NP

MIS LYNX/
C1/V3
SP

MIS LYNX/
C1
WP




EZ-Base System Components

Ratchet key for angled screw channel



MT-ELR10 - Long
MT-ESR10 - Short

Motor key for angled screw channel




MT-ELM10 - Long
MT-ESM10 - Short

EZ-Base prosthetic screw



MN-E0160 - NP
MD-E0220 - SP



A unique screw and driver tip design enable extreme angulation and a wide insertion path.

Digital Restoration

Ti-Base

Anti-rotation

Tightening torque: 30 Ncm

PLATFORM




MIS LYNX/
C1
NP

V3
NP





MIS LYNX/
C1/V3
SP

MIS LYNX/
C1
WP



Ti-Base, 0.5mm



Ti-Base, 1.5mm



Ti-Base, 3mm



Digital Restoration

Ti-Base

Free rotation

Tightening torque: 30 Ncm

i For conical connection, multiple unit, screw-retained restoration, it is recommended to cement the bridge intra-orally.

PLATFORM

Ti-Base, 0.5mm

Ti-Base, 1.5mm

Ti-Base, 3mm




MIS LYNX/
C1/V3

NP






MIS LYNX/
C1/V3

SP



MIS LYNX/
C1

WP



Digital Restoration

Incisor Esthetic Ti Abutments

Anti-rotation

Tightening torque: 30 Ncm

PLATFORM

Incisor esthetic Ti abutment, 0.5mm









Incisor esthetic Ti abutment, 1.5mm

Incisor esthetic Ti abutment, 3mm

MIS LYNX/
C1
NP

V3
NP

MIS LYNX/
C1/V3
SP



Digital Restoration

Digital Model Analog

Tightening torque: 15 Ncm

MIS LYNX/
C1

NP

V3

NP

MIS LYNX/
C1/V3

SP

MIS LYNX/
C1

WP

MULTI-
UNIT

CONNECT
4mm

CONNECT
5.7mm

Ø3.30mm



CN-MAN10

Ø3.30mm



VN-MAN10

Ø3.80mm



CS-MAN10

Ø4.60mm



CW-MAN10

Ø4.80mm



MU-MAN10

Ø4.60mm



MM-MAN40

Ø5.7mm



MM-MAN57

Digital Restoration

Scan Posts

Tightening torque: 15 Ncm

MIS LYNX/
C1
NP

V3
NP

MIS LYNX/
C1/V3
SP

MIS LYNX/
C1
WP

MULTI-
UNIT

CONNECT
4mm

CONNECT
5.7mm

i MIS scan posts are also sold in a kit. CONNECT scan posts are not included in the kit.



MK-0057
SCAN POST KIT,
CONICAL CONNECTION



Digital Restoration
Titanium Blanks
 Anti-rotation

Tightening torque: 30 Ncm

PLATFORM

Titanium blank, Ø9



CN-IB091



VN-IB091



CS-IB091



CW-IB091

Titanium blank, Ø12



CN-IB121



VN-IB121



CS-IB121



CW-IB121

Titanium blank for
 Amann Girrbach, Ø12



CN-AGT12



VN-AGT12



CS-AGT12




CW-AGT12

Milling analogs




CN-MAB10



VN-MAB10



CS-MAB10



CW-MAB10


Cemented Restoration

Cementable Abutments


Tightening torque: 30 Ncm

PLATFORM.


Cementable abutment, 0°



Cementable angulated abutment, 10°



Cementable angulated abutment, 20°



MIS LYNX/
C1
NP

V3
NP

Cemented Restoration

Cementable Abutments


Tightening torque: 30 Ncm

PLATFORM...


MIS LYNX/
C1/V3

SP


Cementable abutment, 0°



Cementable angulated abutment, 15°



Cementable angulated abutment, 25°



Cemented Restoration


Cementable Abutments

Tightening torque: 30 Ncm


PLATFORM..

MIS LYNX/
C1
WP


Cementable abutment, 0°



Cementable angulated abutment, 15°



Cementable angulated abutment, 25°



Cemented Restoration

Cementing posts

Tightening torque: 30 Ncm

PLATFORM.

Cementing post, 1mm

Cementing post, 3mm











MIS LYNX/
C1
NP

V3
NP

MIS LYNX/
C1/V3
SP

MIS LYNX/
C1/V3
SP

MIS LYNX/
C1
WP



Cemented Restoration

CPK - Transgingival Abutments

4mm crown

Tightening torque: 30 Ncm

PLATFORM

Transgingival abutment, 1mm

Transgingival abutment, 2mm


Transgingival abutment, 3mm

Transgingival abutment, 4mm

MIS LYNX/
C1
NP

V3
NP


MIS LYNX/
C1/V3
SP




CN-CPK41




VN-P4140




CS-CPK41



CN-CPK42




VN-P4240




CS-CPK42




CN-CPK43




VN-P4340




CS-CPK43



CN-CPK44



VN-P4440



CS-CPK44

Cemented Restoration

CPK - Transgingival Abutments

6mm crown

Tightening torque: 30 Ncm

PLATFORM.

Transgingival abutment, 1mm

Transgingival abutment, 2mm


Transgingival abutment, 3mm

Transgingival abutment, 4mm


MIS LYNX/
C1
NP




CN-CPK61



CN-CPK62



CN-CPK63




CN-CPK64


V3
NP




VN-P6140



VN-P6240




VN-P6340



VN-P6440


MIS LYNX/
C1/V3
SP




CS-CPK61



CS-CPK62



CS-CPK63



CS-CPK64


MIS LYNX/
C1/V3
SP



VS-P6155



VS-P6255



VS-P6355




VS-P6455


MIS LYNX/
C1
WP




CW-CPK61



CW-CPK62



CW-CPK63



CW-CPK64

Cemented Restoration

CPK -
Transgingival
Abutments

8mm crown

Tightening torque: 30 Ncm

PLATFORM.

Transgingival abutment, 1mm

Transgingival abutment, 2mm

Transgingival abutment, 3mm

Transgingival abutment, 4mm

MIS LYNX/
C1/V3


SP




CS-CPK81



CS-CPK82




CS-CPK83



CS-CPK84

MIS LYNX/
C1/V3

SP




VS-P8155



VS-P8255



VS-P8355



VS-P8455

Cemented Restoration

CPK - Kit Contents



Burnout plastic cap,
anti-rotation




MN-ICH40
MD-ICH40
MW-ICH40

Burnout plastic cap,
free rotation



MN-IC040
MD-IC040
MW-IC040

Plastic healing cap



MM-CHC40
MM-CHC60
MM-CHC80

Transfer plastic cap



MN-IC600
MD-IC800
MW-IC800

Analogs



MN-RSM60
MD-RSM60
MW-RSM60

CPK abutments are also sold in kits

CK-NPK41 ● Ø4, H=1, C=4
VK-NPK41 ● Ø4, H=1, C=4
CK-CPK41 ● Ø4.8, H=1, C=4

CK-NPK42 ● Ø4, H=2, C=4
VK-NPK42 ● Ø4, H=2, C=4
CK-CPK42 ● Ø4.8, H=2, C=4

CK-NPK43 ● Ø4, H=3, C=4
VK-NPK43 ● Ø4, H=3, C=4
CK-CPK43 ● Ø4.8, H=3, C=4

CK-NPK44 ● Ø4, H=4, C=4
VK-NPK44 ● Ø4, H=4, C=4
CK-CPK44 ● Ø4.8, H=4, C=4

CK-NPK61 ● Ø4, H=1, C=6
VK-NPK61 ● Ø4, H=1, C=6
CK-CPK61 ● Ø4.8, H=1, C=6
VK-CPK61 ● Ø4.8, H=1, C=6
CK-WPK61 ● Ø5.5, H=1, C=6

CK-NPK62 ● Ø4, H=2, C=6
VK-NPK62 ● Ø4, H=2, C=6
CK-CPK62 ● Ø4.8, H=2, C=6
VK-CPK62 ● Ø4.8, H=2, C=6
CK-WPK62 ● Ø5.5, H=2, C=6

CK-NPK63 ● Ø4, H=3, C=6
VK-NPK63 ● Ø4, H=3, C=6
CK-CPK63 ● Ø4.8, H=3, C=6
VK-CPK63 ● Ø4.8, H=3, C=6
CK-WPK63 ● Ø5.5, H=3, C=6

CK-NPK64 ● Ø4, H=4, C=6
VK-NPK64 ● Ø4, H=4, C=6
CK-CPK64 ● Ø4.8, H=4, C=6
VK-CPK64 ● Ø4.8, H=4, C=6
CK-WPK64 ● Ø5.5, H=4, C=6

CK-CPK81 ● Ø4.8, H=1, C=8
VK-CPK81 ● Ø5.5, H=1, C=8

CK-CPK82 ● Ø4.8, H=2, C=8
VK-CPK82 ● Ø5.5, H=2, C=8

CK-CPK83 ● Ø4.8, H=3, C=8
VK-CPK83 ● Ø5.5, H=3, C=8
CK-CPK84 ● Ø4.8, H=4, C=8
VK-CPK84 ● Ø5.5, H=4, C=8

Overdenture LOCKiT Attachment System

Tightening torque: 30 Ncm

i All LOCKiT products shown in this catalog are sold subject to individual regulatory certification only.

PLATFORM:

LOCKiT abutment, 1mm

LOCKiT abutment, 2mm


LOCKiT abutment, 3mm

LOCKiT abutment, 4mm

LOCKiT abutment, 5mm

MIS LYNX/
C1/V3


NP



CN-LOC01



CN-LOC02



CN-LOC03



CN-LOC04



CN-LOC05

MIS LYNX/
C1/V3

SP



CS-LOC01



CS-LOC02



CS-LOC03



CS-LOC04



CS-LOC05

MIS LYNX/
C1

WP



CW-LOC01



CW-LOC03



CW-LOC05

Overdenture
LOCKiT
System
Components



LOCKiT abutments
are also sold in kits

- CK-NLC1 ● Ø3.75, H=1
- CK-SLC1 ● Ø3.75, H=1
- CK-WLC1 ● Ø4.5, H=1
- CK-NLC2 ● Ø3.75, H=2
- CK-SLC2 ● Ø3.75, H=2
- CK-NLC3 ● Ø3.75, H=3
- CK-SLC3 ● Ø3.75, H=3
- CK-WLC3 ● Ø4.5, H=3
- CK-NLC4 ● Ø3.75, H=4
- CK-SLC4 ● Ø3.75, H=4
- CK-NLC5 ● Ø3.75, H=5
- CK-SLC5 ● Ø3.75, H=5
- CK-WLC5 ● Ø4.5, H=5

Core tool



MT-LOC10

Ratchet driver



MT-RSD19 - Short
MT-RSD24 - Long

Caps range



MM-BLA10
MM-RED10
MM-GRE10
MM-BLU10
MM-PIN10
MM-CLE10

Male processing
package



MM-LOC02

Metal housing



MM-TIT10

Impression coping



MM-LOC1

Analog



MM-LOA40 - NP/SP
MM-LOA50 - WP

Abutment holder sleeve




MT-LOC11

Overdenture OT-Equators

Tightening torque: 30 Ncm

PLATFORM.

OT-Equator, 1mm




OT-Equator, 2mm



OT-Equator, 3mm



OT-Equator, 4mm



OT-Equator, 5mm



Overdenture

OT-Equators Components

OT-Equators are sold in kits




Impression coping cap



OE-IC002

Analog



OE-RSM60

Metal housing and
transfer caps



OE-MH001
OE-RCB01
OE-RCY01
OE-RCP01
OE-RCW01
OE-RCV01

Handling tool



ET-IT002

Driver for 0.05 inch hex.



MT-RDS30 - Short
MT-RDL30 - Long

Overdenture

Ball Attachment System

Tightening torque: 30 Ncm

PLATFORM.

Ball attachment, 1mm

Ball attachment, 2mm

Ball attachment, 3mm

Ball attachment, 4mm

Ball attachment, 5mm

MIS LYNX/
C1/V3

NP



CN-B1330



CN-B2330




CN-B3330




CN-B4330

MIS LYNX/
C1/V3

SP



CS-B1375



CS-B2375



CS-B3375



CS-B4375




CS-B5375

MIS LYNX/
C1

WP



CW-B1500



CW-B3500

Overdenture

Ball Attachment System Components

Analog




MB-RS375

Titanium ball cap




MB-TTP10

Plastic mounting sleeve



MB-PR010

Gold ball cap



MB-GPS10

Plastic ball cap



MB-SF200

Cap spring



MB-TBS10

Cap spring



MB-TBS20

Standard plastic cap
for metal housing




MB-PPC10

Soft plastic cap for
metal housing




MB-PPR10

Plastic disc



MB-DB235

Double metal housing



MB-DMH10

Ball attachment
anchor key



MT-RB225

Prosthetic Screws and Drivers

Tightening torque: 20/30 Ncm

Prosthetic screw



MN-S0160 - NP
MD-S0220 - SP

Medium prosthetic
screw



MD-S0222

Long prosthetic screw



MD-S0224

Gold prosthetic screw




MD-G0220

CONNECT
prosthetic screw



MM-S0160

Extra-short driver for
0.05 inch hex.



MT-RDM30

Short driver for
0.05 inch hex.



MT-RDS30

Long driver for
0.05 inch hex.




MT-RDL30

Extra-long driver for
0.05 inch hex.




MT-RDE30

Short motor adapter for
0.05 inch hex.



MT-SM005

Long motor adapter for
0.05 inch hex.



MT-LM005

Prosthetic Screws and Keys for Angled Screw Channel

i For use with CONNECT
final esthetic abutment only
and requires angled screw
channel key

Prosthetic screw for
angled screw channel*



MM-SA160

EZ-Base prosthetic screw



MN-E0160 - NP
MD-E0220 - SP




A unique screw and driver tip
design enable extreme angulation
and a wide insertion path.

Ratchet key for angled
screw channel



MT-ELR10 - Long
MT-ESR10 - Short

Motor key for angled
screw channel




MT-ELM10 - Long
MT-ESM10 - Short

MK-0039

Prosthetic Tools Kit

MK-0039
PROSTHETIC TOOLS KIT





All rights reserved. No part of this publication may be reproduced, transcribed, stored in an electronic retrieval system, translated into any language or computer language, or be transmitted in any form whatsoever, without the prior written consent of the publisher. Warning: MIS products referred to in this document should be used by licensed dentists only.

Products may not be available in some regulated markets.
Please contact your local Dentsply Sirona sales office for current product assortment and availability.

